

**Appendix G**

**ASTM E1598**

**Plant Toxicity Test**



## Appendix G

### ASTM E1598 Plant Toxicity Test

#### SOUTHWEST RESEARCH INSTITUTE®

6220 CULEBRA RD. 78208-5198 • P.O. DRAWER 26610 78228-0510 • SAN ANTONIO, TEXAS, USA • (210) 522-5111 • WWW.SWRI.ORG

Chemistry and Chemical Engineering Division  
Department of Analytical and Environmental Chemistry

December 21, 2004

INEEL  
2525 North Fremont  
Idaho Falls, ID 83415-3960


Attn: Kathleen A. Otter

Subject:	INEEL Contract No.:	50157
	INEEL SOW No.:	ER-SOW-473
	SDG Number:	ECR111019A
	SwRI Project Number:	06355.38.00X
	SwRI Task Order Number:	040825-11, 040826-13, 040827-11, 040915-3
	SwRI Sample Receipt Number:	26410, 26431, 26438, 26520
	Samples Received:	August 25, 26, 27, and September 10, 2004
	Plant Toxicity Test	


Dear Ms. Otter:

Enclosed please find the cover page and analytical data for the above referenced project. If you should have any questions, please do not hesitate to call me at (210) 522-5428.

Sincerely,

  
Michael J. Dammann  
Manager

APPROVED:

  
Reza Karimil, Ph.D.  
Director

MJD:mar



DETROIT, MICHIGAN (248) 593-2550 • HOUSTON, TEXAS (713) 577-1577 • WASHINGTON, DC (301) 861-0226





# ASTM E1598

## Plant Toxicity Test

### **Analytical Report**

*Southwest Research Institute*

*Project 06355.38.006*

*Task Orders 040825-11, 040826-13, 040827-11, 040915-3*

### **Prepared by**

*Michael Dammann*

*Radonna Spies*

*Terence O'Brien*

### **Prepared for**

*IDAHO NATIONAL ENGINEERING AND*

*ENVIRONMENTAL LABORATORY*

*ER-SOW-444*

**Sample Identification**

Thirty-four soil samples were received at Southwest Research Institute. The samples were logged in to the Laboratory Information Management System (LIMS), and each sample was assigned a unique identification number, referred to as the SwRI Lab System ID. See Table 1 for cross-reference to the client ID.

**Table 1**

SwRI Lab System ID	Client ID	Date Received
250022	ECR111019A	08/25/04
250023	ECR111029A	08/25/04
250024	ECR112019A	08/25/04
250025	ECR113019A	08/25/04
250026	ECR114019A	08/25/04
250027	ECR115019A	08/25/04
250028	ECR116019A	08/25/04
250029	ECR117019A	08/25/04
250030	ECR118019A	08/25/04
250031	ECR119019A	08/25/04
250032	ECR120019A	08/25/04
250152	ECX111019A	08/26/04
250153	ECX111029A	08/26/04
250154	ECX112019A	08/26/04
250155	ECX113019A	08/26/04
250156	ECX114019A	08/26/04
250157	ECX115019A	08/26/04
250158	ECX116019A	08/26/04
250159	ECX117019A	08/26/04
250160	ECX118019A	08/26/04
250161	ECX119019A	08/26/04
250162	ECX120019A	08/26/04

***Table 1 (cont.)***

SwRI Lab System ID	Client ID	Date Received
250231	ECT111019A	08/27/04
250232	ECT111029A	08/27/04
250233	ECT112019A	08/27/04
250234	ECT113019A	08/27/04
250235	ECT114019A	08/27/04
250236	ECT115019A	08/27/04
250237	ECT116019A	08/27/04
250238	ECT117019A	08/27/04
250239	ECT118019A	08/27/04
250240	ECT119019A	08/27/04
250241	ECT120019A	08/27/04
251336	ECT126019A	09/10/04

### ***Introduction***

For the INEEL Long-Term Ecological Monitoring FY-04 project, the toxicity of site soils to representative seedlings was evaluated using ASTM E1598 *Standard Practice for Determining Early Seedling Growth Tests*. At the conclusion of the 30-day test period, the survival and growth of the test seedlings were compared to a control soil; the results were analyzed statistically for adverse effects.

### ***Test Procedure***

#### ***Control Soils***


- 1) Originally Controls were prepared with – Potting Soil – Packaged by D.C. Organics, Inc. 3101 FM 2331, Godley, TX 76044
- 2) Due to poor performance of the plants in the original control soil, a second set of control samples were started using – Garden Plus® - Enriched Top Soil, Packaged by the Letco Group, Dallas, Texas 15220.

#### ***Test Soils***

Each soil, including the controls, were thoroughly mixed and passed through a 2.0 mm sieve until a sufficient amount was attained to fill the test containers.

Test Trays

Each tray held three 6-pack test container inserts, ten individual containers and two 6-pack empty containers. This was an optimum design since each tray only held one sample with five test species (see Figure 1 below, and Image 2 in Appendix A). This allowed for the use of larger "pots" for better root growth and reduction of possible cross contamination. Thirty-four samples and three duplicates and three controls resulted in a total of forty trays.



**Figure 1. Tray Configuration**

Test Species

The certified organic seeds were purchased from High Mowing Seed Company in Northern Vermont. The company provides lot numbers and germination rates with their seeds. Appendix B includes photocopies of the seed packets, and Appendix A contains images of the actual seeds. Table 2 indicates the five species of plants selected for testing.

**Table 2. Five Plant Species Used for Testing**

Group	Family	Species	Common Name	Lot #	Germ. Rate
Dicotyledons	Cruciferae	Raphanus sativum	Radish	TER8-70	86%
Dicotyledons	Luguminosae	Glycine max	Soybean	MHM9-04	96%
Dicotyledons	Umbelliferae	Daucus carota	Carrot	TER8-03	89%
Monocotyledons	Amaryllidaceae	Allium cepa	Onion	TER9-19	89%
Monocotyledons	Gramineae	Zea mays	Corn	BF8-06.1	94%

*Seed Planting*

The onion, carrot and radish seeds were small; five seeds were planted 1.0–1.5 cm deep in each container. Due to their larger size, three soybean and corn seeds were planted 2.5–4.0 cm deep in each test container. Images 3a-e in Appendix A show the relative size of the seeds. Table 3 is the breakdown of the total number of seeds planted.

**Table 3. Number of Seeds Planted per Replicate**

Plant	Replicates	Seeds/Replicate	Total Seeds
Radish	6	5	30
Soybean	6	3	18
Carrot	6	5	30
Onion	6	5	30
Corn	6	3	18

*Facility*

An indoor growth chamber was prepared in a separate building from the analytical laboratories to prevent toxic contamination or vapors. A table, measuring 2.89m by 2.12m, was made to support the sample planting trays; it allowed space for 40 trays.

*Planting*

The amount of seeds to be planted and then plants to be harvested at the end of the growing period required that the samples and controls be planted in three groupings. One control and one duplicated sample were planted with each grouping of samples. The original control sample planted with each of the planting sets developed a hard crusty layer on the surface that suffocated the plant seeds; this was possibly due to manure present in the soil mix. The initial controls were therefore rejected and new controls were started using a new soil while the samples were still growing. The sample groupings start dates and end dates are listed in Table 11.

*Light Requirements*

The only light source found to produce the minimum light requirements was a 400-Watt metal halide HID lamp. Twelve lamps were positioned in 4 rows of 3 lights over the 40 trays (see Image 1 in Appendix A). The light intensity was monitored daily using a recording photometer, the meter was placed at the distance from the lights equal to that of the soil surface and along the very outside of the growth region. This arrangement leads to the light measurements showing the absolute minimum

exposure the plants would see. The light cycle was controlled using a programmable timer, which turned lights on at 6:30 AM and turned turn off at 8:30 PM. The approximate light intensity averaged 2000fc during the duration of the test. The standard recommends 2000fc, with a range of 300-400  $\mu\text{mol}/\text{m}^2/\text{s}$ . Graph 1 illustrates the light output for the duration of the 30-day test. The height of the lamps was adjusted twice during the 30-days test period to accommodate the growing plants. The initial height adjustment occurred on September 22 with a second height adjustment on September 28. Occasional dips below the 300  $\mu\text{mol}/\text{m}^2/\text{s}$  were due to plant leaves getting between the lights and the photometer.

#### Temperature and Humidity Requirements

The air temperature and humidity were also continuously monitored using a computer program. Graph 2 illustrates the environmental conditions during the 30-day test. The minimum and maximum temperatures for the duration of the test were 18.5 °C and 34.1 °C, respectively. There was one night that the temperature dropped below the ideal minimum temperature of 20°C, due to a cold front. There were multiple occasions that the temperature exceeded the recommended 30°C due to outside temperatures and the heat output from the lights. The humidity ranged from 32.6 to 73.6 %.

#### Watering Schedule

The seeds were initially watered with 20 ml of tap water for the 6-packs and 40 for the individual containers, and then monitored daily. Daily watering was necessary to keep the soil moist. Watering was performed using a hose connected to a metering peristaltic pump. The pump was calibrated to deliver a set amount of water per foot petal trigger. No nutrients were added during the test.

#### Rotating Schedule

To ensure that test species received average uniform lighting, the trays were rotated daily.

- 1) The trays in position A8-A10 and B8-B10 were removed, all remaining trays in Rows A & B were moved to the right three positions, and the trays from positions A8-A10 were placed in positions B1-B3. The trays from B8-B10 were placed in positions A1-A3.
- 2) The trays in position C8-C10 and D8-D10 were removed, all remaining trays in Row C & D were moved to the left three positions, and the trays from positions C8-C10 were placed in positions D1-D3. The trays from D8-D10 were placed in positions C1-C3.

Window	Position 1	Position 2	Position 3	Position 4	Position 5	Position 6	Position 7	Position 8	Position 9	Position 10	ROW D ROW C ROW B ROW A
Window											
Window											

Seedlings

The total number of seedlings emerging for each species per sample was recorded daily. The germination rates are calculated in Table 5a-e.

End of Test

At the conclusion of the study, the above ground plants were harvested and measured. The remaining below ground roots were washed to remove the soil and measured. The above ground seedling and the below ground roots were dried in an oven at  $70 \pm 5^\circ\text{C}$  for a minimum of 48 hours. The number of plants recovered in the harvesting process can be found in Table 5a-e.

Plant Deaths


Table 4 contains the number of the plants that died during the 30-day test. Statistical Analysis Data acquired at the end of the test include seedling emergence data by species for each container, see Tables 6a-e; dry above ground plant weight by species and container, see Tables 7a-e; dry root weight for each species and container, see Table 8a-e; plant height by species and container, see Table 9a-e; and root length by species and container, see Table 10a-e. Each of the above-mentioned tables contains the data for each parameter by container (discrete test replicate) plus the total sum of each of the containers. Each table also contains a sample labeled "Combined" that is the three control sets treated as a single sample. The mean and standard deviation for the combined control sample was used for all statistical

analysis except the calculation of the pooled standard deviation, which used all samples except the combined control. The pooled standard deviation (shown at the bottom of each table) was calculated as the square root of the sum of the squares of deviations from the means divided by the degree of freedom. This pooled STD was used in calculating any data during the T-test analysis. The right hand side of the tables contains statistical information concerning each sample. Number of measurements (n), mathematical average (mean) and standard deviation data were derived in the normal manner on the data from each sample using each container as the individual replicates for calculations. The column labeled "Difference from Control" contains the variance of the sample mean from the combined control mean. The column labeled "Least Significant Difference @ 95%" contains the calculated least significant difference at the single tailed 95% confidence limit using the T-test, based upon the pooled standard deviation for the set, for each sample and the combined control. The final column shows a "Yes" if the sample mean and the control mean are statistically different.

Observations

- 1) A number of plants died during the test period; see Table 4. The some of the plants that died showed signs of wilting around the middle of the test period. Some of the smaller plants (onions, carrots) were soon indistinguishable from the soil and were not harvestable. Most of the rest of the plants had developed to the point were they were recoverable at harvest time and are included in the height, length and weight data.
- 2) Most of the soil samples developed a green tinge (fungus) on the surface of some of the individual containers. This is noted in Table 4.
- 3) In some instances in Table 5a-e, there were discrepancies in the number of plants emerged and number of plants recovered. Some of these differences were attributed to plant deaths; see Table 4. The remaining differences were probably due to loss of plants during the harvesting process. The number of inconsistencies was greatest in the carrots and onions, which were the smallest plants, and the most difficult to separate from the soil. Images 6a and 6b (Appendix A) show the relative size of the radishes, carrots and onions.


Graph 1. Light Out During the 30-day Test



G-13

010008

*Graph 2. Environmental Conditions During the 30-day Test*















010016

Southwest Research Institute

December 16, 2004

Table 6a. CORN Plant Emergence Numbers

Sample ID	Chen ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Standard Deviation	Difference from Control	Least Significant Statistical Difference
250022	ECR111019A	3	3	3	3	2	17	6	2.83	0.41	0.50	0.48 YES
250022D	ECR111019AD	3	2	3	4	3	3	18	6	3.00	0.63	0.67 0.48 YES
250023	ECR111029A	2	2	2	3	3	2	14	6	2.33	0.52	0.00 0.48
250024	ECR112019A	3	3	3	3	3	3	18	6	3.00	0.00	0.67 0.48 YES
250025	ECR113019A	2	3	2	2	3	3	15	6	2.50	0.55	0.17 0.48
250026	ECR114019A	2	3	3	2	3	3	16	6	2.67	0.52	0.33 0.48
250027	ECR115019A	3	3	3	3	3	3	18	6	3.00	0.00	0.67 0.48 YES
250028	ECR116019A	2	3	3	3	3	3	17	6	2.83	0.41	0.50 0.48 YES
250029	ECR117019A	3	3	3	3	2	2	16	6	2.67	0.52	0.33 0.48
250030	ECR118019A	3	3	3	3	3	2	17	6	2.83	0.41	0.50 0.48 YES
250031	ECR119019A	3	3	3	3	3	2	17	6	2.83	0.41	0.50 0.48 YES
250032	ECR120019A	3	3	2	3	3	2	16	6	2.67	0.52	0.33 0.48
250152	ECX111019A	3	2	3	2	3	3	16	6	2.67	0.52	0.33 0.48
250152D	ECX111019AD	3	3	3	3	3	2	17	6	2.83	0.41	0.50 0.48 YES
250153	ECX111029A	3	3	3	3	3	3	18	6	3.00	0.00	0.67 0.48 YES
250154	ECX112019A	3	2	3	3	3	3	17	6	2.83	0.41	0.50 0.48 YES
250155	ECX113019A	3	3	2	2	3	3	16	6	2.67	0.52	0.33 0.48
250156	ECX114019A	3	3	2	3	3	3	17	6	2.83	0.41	0.50 0.48 YES
250157	ECX115019A	2	2	3	2	2	3	14	6	2.33	0.52	0.00 0.48
250158	ECX116019A	3	3	3	3	2	2	16	6	2.67	0.52	0.33 0.48
250159	ECX117019A	3	3	3	3	3	3	18	6	3.00	0.00	0.67 0.48 YES
250160	ECX118019A	3	2	3	3	3	3	17	6	2.83	0.41	0.50 0.48 YES
250161	ECX119019A	3	3	3	3	3	2	17	6	2.83	0.41	0.50 0.48 YES
250162	ECX120019A	3	3	3	1	2	2	14	6	2.33	0.82	0.00 0.48
250231	ECT111019A	3	3	3	3	3	3	18	6	3.00	0.00	0.67 0.48 YES
250232	ECT111029A	2	2	3	3	3	3	16	6	2.67	0.32	0.33 0.48
250233	ECT112019A	3	3	1	2	3	3	15	6	2.50	0.84	0.17 0.48
250234	ECT113019A	2	3	3	3	3	3	17	6	2.83	0.41	0.50 0.48 YES
250235	ECT114019A	1	1	1	2	2	3	10	6	1.67	0.82	-0.67 0.48 YES
250236	ECT115019A	0	0	0	0	0	0	0	6	0.00	0.00	-2.33 0.48 YES
250237	ECT116019A	0	3	1	2	1	0	7	6	1.17	1.17	-1.17 0.48 YES
250238	ECT117019A	3	3	0	1	2	2	11	6	1.83	1.17	-0.50 0.48 YES
250239	ECT118019A	1	1	0	0	1	0	3	6	0.50	0.55	-1.83 0.48 YES
250240	ECT119019A	0	0	0	0	0	0	0	6	0.00	0.00	-2.33 0.48 YES
250241	ECT120019A	2	0	3	0	0	0	5	6	0.83	1.33	-1.50 0.48 YES
251336	ECT126019A	2	2	2	2	2	3	13	6	2.17	0.41	-0.17 0.48
251336D	ECT126019AD	3	1	1	3	3	1	12	6	2.00	1.10	-0.33 0.48
Control-1R		3	2	2	2	2	3	14	6	2.33	0.52	0.00 0.48
Control-2R		2	3	3	2	3	2	15	6	2.50	0.55	0.17 0.48
Control-3R		2	2	3	3	2	1	13	6	2.17	0.75	-0.17 0.48
Combined								42	18	2.33	0.59	0.00 0.33
Pooled STD		0.59										



010018

*Southwest Research Institute*

*December 16, 2004*

*Table 6c. CARROT Plant Emergence Numbers*

Sample ID	Chen ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR11019A	5	5	5	6	4	4	29	6	4.83	0.75	2.94	0.78	YES
250022D	ECR11019AD	4	5	5	4	6	5	29	6	4.83	0.75	2.94	0.78	YES
250023	ECR11029A	4	5	4	4	4	5	26	6	4.33	0.52	2.44	0.78	YES
250024	ECR112019A	5	4	4	4	5	5	27	6	4.50	0.55	2.61	0.78	YES
250025	ECR113019A	5	4	5	5	4	4	27	6	4.50	0.55	2.61	0.78	YES
250026	ECR114019A	5	5	5	5	5	4	29	6	4.83	0.41	2.94	0.78	YES
250027	ECR115019A	5	4	5	5	4	4	27	6	4.50	0.55	2.61	0.78	YES
250028	ECR116019A	1	5	5	5	5	5	26	6	4.33	1.63	2.44	0.78	YES
250029	ECR117019A	3	2	3	3	4	5	20	6	3.33	1.03	1.44	0.78	YES
250030	ECR118019A	4	4	3	4	5	5	25	6	4.17	0.75	2.28	0.78	YES
250031	ECR119019A	1	3	3	2	4	3	16	6	2.67	1.03	0.78	0.78	
250032	ECR120019A	0	0	0	0	0	0	0	6	0.00	0.00	-1.89	0.78	YES
250152	ECX111019A	0	4	3	2	5	1	15	6	2.50	1.87	0.61	0.78	
250152D	ECX111019AD	5	5	4	5	3	5	27	6	4.50	0.84	2.61	0.78	YES
250153	ECX111029A	4	4	5	5	5	4	27	6	4.50	0.55	2.61	0.78	YES
250154	ECX112019A	4	0	4	4	0	6	18	6	3.00	2.45	1.11	0.78	YES
250155	ECX113019A	2	4	4	4	5	2	21	6	3.50	1.22	1.61	0.78	YES
250156	ECX114019A	5	4	5	5	5	5	29	6	4.83	0.41	2.94	0.78	YES
250157	ECX115019A	0	1	1	0	1	1	4	6	0.67	0.52	-1.22	0.78	YES
250158	ECX116019A	4	4	3	2	3	5	21	6	3.50	1.05	1.61	0.78	YES
250159	ECX117019A	4	5	5	3	3	3	23	6	3.83	0.98	1.94	0.78	YES
250160	ECX118019A	5	4	4	4	5	4	26	6	4.33	0.52	2.44	0.78	YES
250161	ECX119019A	4	2	5	5	4	4	24	6	4.00	1.10	2.11	0.78	YES
250162	ECX120019A	5	5	2	5	4	5	26	6	4.33	1.21	2.44	0.78	YES
250231	ECT111019A	4	2	4	0	5	4	19	6	3.17	1.83	1.28	0.78	YES
250232	ECT111029A	0	3	4	2	1	1	11	6	1.83	1.47	-0.06	0.78	
250233	ECT112019A	5	5	3	4	4	4	25	6	4.17	0.75	2.28	0.78	YES
250234	ECT113019A	4	4	5	5	4	5	27	6	4.50	0.55	2.61	0.78	YES
250235	ECT114019A	0	1	0	0	1	0	2	6	0.33	0.52	-1.56	0.78	YES
250236	ECT115019A	2	0	0	0	0	0	2	6	0.33	0.62	-1.56	0.78	YES
250237	ECT116019A	0	0	0	0	0	0	0	6	0.00	0.00	-1.89	0.78	YES
250238	ECT117019A	0	0	0	0	0	1	1	6	0.17	0.41	-1.72	0.78	YES
250239	ECT118019A	0	0	0	0	0	1	0	6	0.17	0.41	-1.72	0.78	YES
250240	ECT119019A	0	0	0	0	0	0	0	6	0.00	0.00	-1.89	0.78	YES
250241	ECT120019A	0	0	0	0	0	1	1	6	0.17	0.41	-1.72	0.78	YES
251336	ECT126019A	0	0	0	1	0	0	1	6	0.17	0.41	-1.72	0.78	YES
251336D	ECT126019AD	0	0	1	1	0	0	2	6	0.33	0.52	-1.56	0.78	YES
Control-1R		4	3	3	1	2	0	13	6	2.17	1.47	0.28	0.78	
Control-2R		5	3	3	3	2	2	18	6	3.00	1.10	1.11	0.78	
Control-3R		1	0	1	0	0	1	3	6	0.50	0.55	-1.39	0.78	
Combined								34	18	1.89	1.49	0.00	0.54	
Pooled STD		0.96												

010019

Southwest Research Institute

December 16, 2004

*Table 6d. RADISH Plant Emergence Numbers*

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Last Significant Difference @ 95% Statistical Significance
250022	ECR111019A	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67	
25002D	ECR111019AD	5	5	5	5	4	29	6	4.83	0.41	0.06	0.67	
250023	ECR111029A	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67	
250024	ECR112019A	5	5	6	5	5	31	6	5.17	0.41	0.39	0.67	
250025	ECR113019A	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67	
250026	ECR114019A	5	5	4	5	5	29	6	4.83	0.41	0.06	0.67	
250027	ECR115019A	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67	
250028	ECR116019A	4	5	6	5	4	5	29	6	4.83	0.75	0.06	0.67
250029	ECR117019A	5	4	5	5	4	6	29	6	4.83	0.75	0.06	0.67
250030	ECR118019A	4	5	4	3	5	5	26	6	4.33	0.62	-0.44	0.67
250031	ECR119019A	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67	
250032	ECR120019A	3	5	3	2	5	2	20	6	3.63	1.37	-1.44	0.67
250152	ECX111019A	1	4	0	2	5	4	16	6	2.67	1.97	-2.11	0.67
250152D	ECX111019AD	5	3	5	5	5	30	6	5.00	0.00	0.22	0.67	
250153	ECX111029A	5	5	5	4	5	5	29	6	4.83	0.41	0.06	0.67
250154	ECX112019A	1	5	0	4	4	2	16	6	2.67	1.97	-2.11	0.67
250155	ECX113019A	5	3	5	4	5	4	26	6	4.33	0.82	-0.44	0.67
250156	ECX114019A	5	6	4	4	5	5	29	6	4.83	0.75	0.06	0.67
250157	ECX115019A	0	3	0	0	0	0	3	6	0.50	1.22	-4.28	0.67
250158	ECX116019A	5	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67
250159	ECX117019A	4	5	5	5	6	4	29	6	4.83	0.75	0.06	0.67
250160	ECX118019A	5	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67
250161	ECX119019A	5	5	5	5	5	4	29	6	4.83	0.41	0.06	0.67
250162	ECX120019A	4	5	5	5	4	4	27	6	4.50	0.55	-0.28	0.67
250231	ECT111019A	5	4	5	4	5	4	27	6	4.50	0.55	-0.28	0.67
250232	ECT111029A	5	5	5	5	5	3	28	6	4.67	0.82	-0.11	0.67
250233	ECT112019A	5	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67
250234	ECT113019A	5	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67
250235	ECT114019A	4	5	4	5	5	4	27	6	4.50	0.55	-0.28	0.67
250236	ECT115019A	5	2	4	2	3	3	19	6	3.17	1.17	-1.61	0.67
250237	ECT116019A	4	5	5	2	4	1	21	6	3.50	1.64	-1.28	0.67
250238	ECT117019A	5	4	4	3	5	5	26	6	4.33	0.82	-0.44	0.67
250239	ECT118019A	4	4	0	2	4	1	15	6	2.50	1.76	-2.28	0.67
250240	ECT119019A	1	0	1	1	1	1	5	6	0.83	0.41	-3.94	0.67
250241	ECT120019A	4	2	3	4	3	1	14	6	2.33	1.21	-2.44	0.67
251336	ECT126019A	5	5	5	4	5	4	28	6	4.67	0.52	-0.11	0.67
251336D	ECT126019AD	5	5	5	5	5	4	29	6	4.83	0.41	0.06	0.67
Control-1R		5	4	5	5	5	5	29	6	4.83	0.41	0.06	0.67
Control-2R		5	5	5	5	5	5	30	6	5.00	0.00	0.22	0.67
Control-3R		5	5	4	4	4	5	27	6	4.50	0.55	-0.28	0.67
Combined								86	18	4.78	0.43	0.00	0.47

Pooled STD 0.83

**010020***Southwest Research Institute*

December 16, 2004

**Table 6e. SOYBEAN Plant Emergence Numbers**

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Standard Deviation	Difference from Control	Lead Significance @ 95%	Statistical Difference	
250022	ECR11019A	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51		
250022D	ECR11019AD	2	3	3	3	3	17	6	2.83	0.41	0.28	0.51		
250023	ECR11029A	3	3	3	3	0	15	6	2.50	1.22	-0.06	0.51		
250024	ECR112019A	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51		
250025	ECR113019A	3	0	3	3	2	3	14	6	2.43	1.21	-0.22	0.51	
250026	ECR114019A	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51		
250027	ECR115019A	3	3	3	3	2	3	17	6	2.83	0.41	0.28	0.51	
250028	ECR116019A	2	3	3	3	3	3	17	6	2.83	0.41	0.28	0.51	
250029	ECR117019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250030	ECR118019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250031	ECR119019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250032	ECR120019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250152	ECK111019A	3	3	3	3	3	2	17	6	2.83	0.41	0.28	0.51	
250152D	ECX11019AD	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250153	ECK111029A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250154	ECX112019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250155	ECX113019A	1	3	3	3	3	3	16	6	2.67	0.82	0.11	0.51	
250156	ECX114019A	2	3	3	3	3	3	17	6	2.83	0.41	0.28	0.51	
250157	ECX115019A	2	1	1	1	2	3	10	6	1.67	0.82	-0.89	0.51	YES
250158	ECX116019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250159	ECX117019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250160	ECX118019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250161	ECX119019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250162	ECX120019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250231	ECT111019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250232	ECT111029A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250233	ECT112019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250234	ECT113019A	3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
250235	ECT114019A	1	3	2	3	0	3	12	6	2.00	1.26	-0.56	0.51	YES
250236	ECT115019A	1	3	2	3	1	0	10	6	1.67	1.21	-0.89	0.51	YES
250237	ECT116019A	3	0	2	1	0	2	8	6	1.33	1.21	-1.22	0.51	YES
250238	ECT117019A	2	3	3	3	2	3	16	6	2.67	0.52	0.11	0.51	
250239	ECT118019A	0	1	0	0	0	0	1	6	0.17	0.41	-2.39	0.51	YES
250240	ECT119019A	0	0	0	0	0	0	0	6	0.00	0.00	-2.56	0.51	YES
250241	ECT120019A	3	1	2	3	0	0	9	6	1.50	1.38	-1.06	0.51	YES
251336	ECT126019A	3	3	3	3	1	3	16	6	2.67	0.82	0.11	0.51	
251336D	ECT126019AD	3	2	3	3	2	1	14	6	2.33	0.82	-0.22	0.51	
Control-1R		3	3	3	2	3	3	17	6	2.83	0.41	0.28	0.51	
Control-2R		3	3	3	3	3	3	18	6	3.00	0.00	0.44	0.51	
Control-3R		3	3	0	3	2	0	11	6	1.83	1.47	-0.72	0.51	
Combined								46	18	2.56	0.98	0.00	0.35	
Pooled STD		0.63												

**010021**

*Southwest Research Institute*

*December 16, 2004*

*Table 7a. CORN Total Plant Weights (g)*

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	0.9600	0.9434	1.2232	1.3908	0.9672	1.0294	6.5140	6	1.086	0.182	-0.120	0.204 YES
250022D	ECR111019AD	2.0863	1.2393	1.5399	1.0965	1.2389	1.3634	8.5843	6	1.431	0.357	-0.665	0.204 YES
250023	ECR111029A	0.9167	0.9437	0.9388	1.3514	1.3233	1.1268	6.6007	6	1.100	0.199	-0.334	0.204 YES
250024	ECR112019A	1.2369	1.7727	0.9956	1.5537	1.2811	1.2800	8.1200	6	1.353	0.271	-0.587	0.204 YES
250025	ECR113019A	0.5868	1.1794	0.7745	0.6250	1.0343	1.3059	5.5059	6	0.914	0.299	-0.152	0.204
250026	ECR114019A	0.8933	1.6393	1.0949	0.8968	1.2018	1.1755	6.9016	6	1.150	0.274	-0.384	0.204 YES
250027	ECR115019A	1.5166	1.1907	1.1965	1.1079	1.4786	1.5892	8.0795	6	1.347	0.204	-0.581	0.204 YES
250028	ECR116019A	0.7608	1.4029	1.2381	0.8766	1.3155	0.9309	6.5248	6	1.087	0.265	-0.371	0.204 YES
250029	ECR117019A	1.4720	0.7590	1.5667	1.5325	1.1219	1.0148	7.2669	6	1.211	0.299	-0.445	0.204 YES
250030	ECR118019A	0.8683	0.8267	1.2368	1.4256	1.1382	1.0132	6.5088	6	1.095	0.228	-0.319	0.204 YES
250031	ECR119019A	0.9906	0.7623	1.3669	1.3029	1.6136	1.2259	7.2610	6	1.210	0.298	-0.444	0.204 YES
250032	ECR120019A	0.8813	0.9063	0.8768	1.1639	1.0712	0.8736	5.7731	6	0.962	0.124	-0.196	0.204
250152	ECX111019A	0.5925	0.5412	0.9179	0.5540	1.0937	1.0481	4.7474	6	0.791	0.258	-0.025	0.204
250152D	ECX111019AD	1.3134	1.2085	0.7700	1.1619	1.3281	1.0448	6.8267	6	1.138	0.208	-0.372	0.204 YES
250153	ECX111029A	0.8299	1.1366	1.4039	1.2063	0.5155	1.4488	6.5410	6	1.090	0.358	-0.124	0.204 YES
250154	ECX112019A	0.8674	0.8715	1.2130	0.6417	1.1613	0.8503	5.6052	6	0.934	0.215	-0.168	0.204
250155	ECX113019A	1.2731	0.9024	1.0873	1.3829	1.2738	1.6426	7.5621	6	1.260	0.253	-0.494	0.204 YES
250156	ECX114019A	0.7874	0.9861	0.7766	1.2342	1.2550	1.2209	6.2602	6	1.043	0.225	-0.277	0.204 YES
250157	ECX115019A	0.6736	0.8106	1.1047	1.0242	1.1065	1.3292	6.0488	6	1.008	0.234	-0.242	0.204 YES
250158	ECX116019A	0.6406	1.2196	1.0404	1.2236	0.9355	0.9359	6.0016	6	1.000	0.216	-0.234	0.204 YES
250159	ECX117019A	0.9889	1.1366	1.0659	1.3783	1.4587	1.3294	7.2099	6	1.202	0.211	-0.436	0.204 YES
250160	ECX118019A	1.3179	1.0991	1.2635	1.3968	1.6118	1.5709	8.2600	6	1.377	0.193	-0.611	0.204 YES
250161	ECX119019A	1.0308	0.4236	1.0810	1.3715	1.1124	0.9479	5.9673	6	0.995	0.314	-0.229	0.204 YES
250162	ECX120019A	1.2025	1.4898	1.1139	0.8653	0.7367	1.0936	6.5018	6	1.084	0.264	-0.318	0.204 YES
250231	ECT111019A	1.3600	1.4035	1.3660	1.2520	1.1691	1.2346	7.7852	6	1.298	0.992	-0.531	0.204 YES
250232	ECT111029A	0.7854	0.8595	1.0549	0.9968	0.9052	0.8482	5.4500	6	0.908	0.100	-0.142	0.204
250233	ECT112019A	1.1782	1.4113	0.8727	1.2601	0.9388	1.4800	7.1411	6	1.190	0.246	-0.424	0.204 YES
250234	ECT113019A	1.0364	1.1600	1.2901	1.9296	1.6774	1.1820	8.7321	6	1.455	0.340	-0.689	0.204 YES
250235	ECT114019A	0.5105	0.4834	0.2604	0.9664	0.6664	0.6616	2.9487	6	0.491	0.161	-0.275	0.204
250236	ECT115019A							0.0000	6	0.000	0.000	-0.766	N/A YES
250237	ECT116019A		0.7662	0.4870	0.5640	0.1577		1.9749	4	0.494	0.253	-0.272	0.240 YES
250238	ECT117019A	0.9146	1.1415		0.6614	0.4403	0.6810	3.8388	5	0.766	0.268	0.002	0.219
250239	ECT118019A		0.5623			0.3819		0.9442	2	0.472	0.128	-0.294	0.326
250240	ECT119019A							0.0000	6	0.000	0.000	-0.766	N/A YES
250241	ECT120019A	0.2950		1.2330				1.5289	2	0.764	0.664	-0.002	0.326
251316	ECT126019A	0.7658	1.0253	0.8200	1.0727	0.8680	1.0344	5.5862	6	0.931	0.129	0.165	0.204
251336D	ECT126019AD	0.8883	0.5719	0.4699	0.8768	0.9813	0.5103	4.2985	6	0.716	0.223	-0.050	0.204
Control-1R		0.6420	0.5013	0.0736	0.6509	0.8219	0.8574	3.5471	6	0.591	0.285	-0.175	0.204
Control-2R		1.0569	1.2322	0.9848	0.5720	1.2483	0.8347	5.9289	6	0.988	0.256	-0.222	0.204
Control-3R		0.8416	0.8071	1.1973	0.7975	0.3489	0.3202	4.3126	6	0.719	0.333	-0.047	0.204
Combined								13.7886	18	0.766	0.324	0.000	0.142

Pooled STD 0.25

010022

Southwest Research Institute

December 16, 2004

Table 7b. ONION Total Plant Weights (g)

Sample ID	Plant ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	0.0370	0.0433	0.0222	0.0444	0.0389	0.0680	0.2538	6	0.042	0.015	0.029	-0.013	YES
250023D	ECR111019AD	0.0352	0.0444	0.0419	0.0303	0.0764	0.0643	0.2935	6	0.049	0.018	0.035	-0.013	YES
250023	ECR111029A	0.0215		0.0318	0.0223	0.0307	0.0535	0.1598	5	0.032	0.013	0.018	-0.014	YES
250024	ECR112019A	0.0615	0.0787	0.0959	0.0742	0.0404	0.0821	0.4328	6	0.072	0.019	0.059	-0.013	YES
250025	ECR113019A	0.0231	0.0216	0.0287	0.0396	0.0506	0.0298	0.1934	6	0.032	0.011	0.019	-0.013	YES
250026	ECR114019A	0.0304	0.0286	0.0568	0.0400	0.0466	0.0489	0.2513	6	0.047	0.011	0.028	-0.013	YES
250027	ECR115019A	0.0308	0.0613	0.0666	0.0560	0.0544	0.0730	0.3421	6	0.057	0.015	0.043	-0.013	YES
250028	ECR116019A	0.0135	0.0195	0.0181	0.0676	0.0613	0.0156	0.1956	6	0.033	0.025	0.019	-0.013	YES
250029	ECR117019A	0.0473		0.0188	0.0173	0.0211	0.0082	0.1127	5	0.023	0.015	0.009	-0.014	
250030	ECR118019A	0.0271	0.0306	0.0287	0.0595	0.0255	0.0682	0.2396	6	0.049	0.019	0.026	-0.013	YES
250031	ECR119019A	0.0182	0.0374		0.0434	0.0569		0.1559	4	0.039	0.016	0.025	-0.015	YES
250032	ECR120019A	0.0351	0.0084	0.0483	0.0219	0.0889		0.2026	5	0.041	0.011	0.027	-0.014	YES
250152	ECX111019A	0.0218	0.0342	0.0014	0.0020	0.0583		0.1267	5	0.025	0.025	0.012	-0.014	
250152D	ECX111019AD	0.0286	0.0558	0.0297	0.0202	0.0498	0.0122	0.1963	6	0.033	0.017	0.019	-0.013	YES
250153	ECX111029A		0.0085	0.0219		0.0089	0.0186	0.0579	4	0.014	0.007	0.001	-0.015	
250154	ECX112019A	0.0159	0.0039		0.0200	0.0038	0.0200	0.0636	5	0.013	0.008	-0.001	-0.014	
250155	ECX113019A	0.0039	0.0514		0.0048	0.0300	0.0361	0.1262	5	0.025	0.021	0.012	-0.014	
250156	ECX114019A	0.0425	0.0349	0.0368	0.0320	0.0441	0.0309	0.2212	6	0.037	0.005	0.023	-0.013	YES
250157	ECX115019A		0.0099					0.0099	1	0.010	0.000	-0.004	0.028	
250158	ECX116019A	0.0193	0.0401	0.0166	0.0246	0.0498	0.0555	0.2059	6	0.034	0.016	0.021	-0.015	YES
250159	ECX117019A	0.0322	0.0327	0.0526	0.0084	0.0418	0.0139	0.1816	6	0.040	0.017	0.017	-0.013	YES
250160	ECX118019A	0.0323	0.0412	0.0451	0.0523	0.0440	0.0402	0.2551	6	0.043	0.007	0.029	-0.013	YES
250161	ECX119019A	0.0232	0.0205	0.0282	0.0180	0.0073	0.0325	0.1297	6	0.022	0.009	0.008	-0.013	
250162	ECX120019A	0.0388	0.0241	0.0550	0.0446	0.0346	0.0580	0.2551	6	0.043	0.013	0.029	-0.013	YES
250231	ECT111019A		0.0049		0.0169	0.0397	0.0348	0.0963	4	0.024	0.016	0.011	-0.015	
250232	ECT111029A	0.0135		0.0250	0.0225	0.0146		0.0756	4	0.019	0.006	0.003	-0.015	
250233	ECT112019A	0.0195	0.0205	0.0239	0.0204	0.0092	0.0376	0.1311	6	0.022	0.009	0.008	-0.013	
250234	ECT113019A	0.0386	0.0569	0.0134	0.0220	0.0328	0.0268	0.1905	6	0.032	0.015	0.018	-0.013	YES
250235	ECT114019A		0.0265	0.0117	0.0124			0.0506	3	0.017	0.008	0.003	-0.017	
250236	ECT115019A		0.0389		0.0092	0.0032		0.0513	4	0.017	0.019	0.004	-0.017	
250237	ECT116019A							0.0000	0	0.000	0.000	-0.014	NA	YES
250238	ECT117019A	0.0350			0.0819	0.0644	0.0395	0.2208	4	0.055	0.022	0.042	-0.015	YES
250239	ECT118019A					0.0209		0.0209	1	0.021	0.000	0.007	0.028	
250240	ECT119019A			0.0065				0.0065	1	0.007	0.000	-0.007	0.028	
250241	ECT120019A				0.0036	0.0193		0.0229	2	0.011	0.011	-0.002	0.020	
251336	ECT126019A			0.0112				0.0112	1	0.011	0.000	-0.002	0.028	
251336D	ECT126019AD	0.0060		0.0109			0.0103	0.0272	3	0.009	0.003	-0.004	0.017	
Control-1R		0.0140	0.0104	0.0080	0.0082	0.0152	0.0231	0.0789	6	0.013	0.006	0.000	0.013	
Control-2R		0.0086	0.0136	0.0369		0.0130		0.0721	4	0.018	0.013	0.004	0.015	
Control-3R		0.0090	0.0076	0.0081	0.0083		0.0199	0.0523	5	0.010	0.005	-0.003	0.014	
Combined								0.2033	15	0.014	0.008	0.000	0.009	
Pooled STD		0.015												

010023

Southwest Research Institute

December 16, 2004

Table 7c. CARROT Total Plant Weights (g)

Sample ID	Cheli ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95% Statistical Difference
250022	ECR111019A	0.0522	0.0685	0.0667	0.0732	0.0555	0.0503	0.3664	6	0.051	0.019	0.055	0.018 YES
25002D	ECR111019AD	0.0066	0.0528	0.0573	0.0697	0.0799	0.0743	0.3406	6	0.057	0.027	0.051	0.018 YES
250023	ECR111029A	0.0292	0.0299	0.0378	0.0417	0.0288	0.0542	0.2216	6	0.047	0.010	0.031	0.018 YES
250024	ECR112019A	0.0580	0.0505	0.0728	0.0647	0.0734	0.0707	0.3901	6	0.055	0.009	0.059	0.018 YES
250025	ECR113019A	0.0790	0.0368	0.0487	0.0640	0.0499	0.0511	0.3295	6	0.055	0.015	0.049	0.018 YES
250026	ECR114019A	0.0529	0.0288	0.0578	0.0313	0.0397	0.0487	0.2592	6	0.043	0.012	0.037	0.018 YES
250027	ECR115019A	0.0511	0.0618	0.0526	0.0850	0.0647	0.0550	0.3702	6	0.052	0.013	0.056	0.018 YES
250028	ECR116019A	0.0129	0.1172	0.0933	0.0718	0.0920	0.0685	0.4557	6	0.076	0.046	0.070	0.018 YES
250029	ECR117019A	0.0270	0.0231	0.1235	0.0833	0.0677	0.0690	0.3936	6	0.066	0.017	0.060	0.018 YES
250030	ECR118019A	0.0839	0.0799	0.0665	0.0764	0.1080	0.0887	0.5034	6	0.084	0.014	0.078	0.018 YES
250031	ECR119019A	0.0926	0.1255	0.1066	0.0614	0.1288	0.1231	0.6380	6	0.106	0.026	0.100	0.018 YES
250032	ECR120019A							0.0000	0	0.000	0.000	-0.006	NA YES
250152	ECX111019A		0.1034	0.0353	0.0678	0.1292	0.0142	0.3499	5	0.070	0.051	0.054	0.020 YES
250152D	ECX111019AD	0.0738	0.0738	0.0417	0.0601	0.0302	0.0759	0.3555	6	0.059	0.019	0.053	0.018 YES
250153	ECX111029A	0.0626	0.0443	0.1138	0.1003	0.0600	0.0878	0.4688	6	0.078	0.027	0.072	0.018 YES
250154	ECX112019A	0.0349		0.0008	0.0272		0.0140	0.0769	4	0.019	0.015	0.013	0.022
250155	ECX113019A	0.0907	0.1326	0.1416	0.1118	0.1567	0.1062	0.7395	6	0.124	0.025	0.117	0.018 YES
250156	ECX114019A	0.0579	0.0706	0.0778	0.0591	0.0880	0.0356	0.4090	6	0.068	0.013	0.062	0.018 YES
250157	ECX115019A		0.0045	0.0042		0.0180	0.0012	0.0279	4	0.007	0.007	0.001	0.022
250158	ECX116019A	0.0743	0.0727	0.0841	0.0568	0.0469	0.1101	0.4459	6	0.074	0.022	0.068	0.018 YES
250159	ECX117019A	0.0330	0.1009	0.0617	0.0307	0.0528	0.0409	0.3200	6	0.053	0.026	0.047	0.018 YES
250160	ECX118019A	0.0800	0.0622	0.0651	0.0330	0.0751	0.0631	0.3785	6	0.063	0.016	0.057	0.018 YES
250161	ECX119019A	0.0381	0.0304	0.0176	0.0593	0.0570	0.0329	0.2353	6	0.039	0.016	0.033	0.018 YES
250162	ECX120019A	0.0774	0.0751	0.0354	0.0818	0.0690	0.0699	0.4096	6	0.068	0.017	0.062	0.018 YES
250231	ECT111019A	0.0532	0.0784	0.0940		0.0943	0.0889	0.4088	5	0.082	0.037	0.076	0.020 YES
250232	ECT111029A	0.0538	0.0563	0.0418	0.0272	0.0282	0.2073	0.2073	5	0.041	0.021	0.035	0.020 YES
250233	ECT112019A	0.0580	0.0552	0.0308	0.0215	0.0422	0.0305	0.2382	6	0.040	0.015	0.034	0.018 YES
250234	ECT113019A	0.0573	0.0512	0.0630	0.0573	0.0419	0.0731	0.3438	6	0.057	0.011	0.031	0.018 YES
250235	ECT114019A		0.0080			0.0638		0.0718	2	0.036	0.039	0.030	0.029 YES
250236	ECT115019A	0.0152						0.0152	1	0.015	0.006	0.009	0.040
250237	ECT116019A							0.0000	0	0.000	0.000	-0.006	NA YES
250238	ECT117019A							0.0031	1	0.003	0.000	-0.003	0.040
250239	ECT118019A					0.0009		0.0009	1	0.001	0.000	-0.005	0.040
250240	ECT119019A							0.0000	0	0.000	0.000	-0.006	NA YES
250241	ECT120019A							0.0123	1	0.012	0.009	0.006	0.040
251336	ECT126019A				0.0521			0.0521	1	0.052	0.000	0.046	0.040 YES
251336D	ECT126019AD			0.0019	0.0401			0.0420	2	0.021	0.027	0.015	0.029
Control-1R		0.0106	0.0042	0.0065	0.0039	0.0039		0.0291	5	0.006	0.004	0.000	0.020
Control-2R		0.0109	0.0084	0.0053	0.0094	0.0047	0.0083	0.0470	6	0.008	0.002	0.002	0.018
Control-3R		0.0046	0.0022			0.0012	0.0080	3	0.003	0.002	-0.003	0.024	
Combined								0.0841	14	0.006	0.003	0.000	0.014
Pooled STD		0.022											

010024

Southwest Research Institute

December 16, 2004

Table 7d. RADISH Total Plant Weights (g)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Last Significant Difference @ 95%	Statistical difference
250022	ECR111019A	0.2310	0.2446	0.2399	0.2275	0.3184	0.2890	1.5504	6	0.258	0.037	0.067	0.045	YES
25002D	ECR111019AD	0.3210	0.2753	0.3373	0.3009	0.2475	0.2737	1.7557	6	0.293	0.033	0.101	0.045	YES
250023	ECR111029A	0.2567	0.3637	0.3758	0.2383	0.2826	0.3112	1.8283	6	0.305	0.056	0.113	0.045	YES
250024	ECR112019A	0.2942	0.3367	0.3345	0.3066	0.3421	0.2479	1.8620	6	0.310	0.046	0.119	0.045	YES
250025	ECR113019A	0.3213	0.2629	0.2840	0.2737	0.2482	0.2729	1.6630	6	0.277	0.025	0.086	0.045	YES
250026	ECR114019A	0.2061	0.2252	0.2410	0.1827	0.2284	0.2336	1.3170	6	0.220	0.021	0.024	0.045	
250027	ECR115019A	0.2713	0.3030	0.2270	0.2332	0.2534	0.2311	1.5190	6	0.253	0.030	0.062	0.045	YES
250028	ECR116019A	0.1961	0.2929	0.3354	0.2737	0.2119	0.2799	1.5899	6	0.265	0.052	0.074	0.045	YES
250029	ECR117019A	0.4942	0.3649	0.3996	0.3566	0.2908	0.3228	2.2289	6	0.371	0.071	0.180	0.045	YES
250030	ECR118019A	0.2704	0.2931	0.2149	0.2237	0.3292	0.2605	1.5918	6	0.265	0.043	0.024	0.045	YES
250031	ECR119019A	0.4813	0.3110	0.3875	0.3760	0.3119	0.4235	2.2915	6	0.382	0.066	0.190	0.045	YES
250032	ECR120019A	0.2027	0.3110	0.3346	0.1552	0.3023	0.2254	1.5312	6	0.255	0.071	0.064	0.045	YES
250152	ECX111019A	0.2568	0.5939	0.3754	0.4449	0.4904	2.1614	5	0.432	0.126	0.241	0.048	YES	
250152D	ECX111019AD	0.2336	0.2838	0.2553	0.2333	0.2170	0.2587	1.4817	6	0.247	0.024	0.056	0.045	YES
250153	ECX111029A	0.2658	0.3358	0.4687	0.3418	0.2833	0.3411	2.0365	6	0.359	0.071	0.148	0.045	YES
250154	ECX112019A	0.1957	0.2875	0.2875	0.1530	0.6362	1	0.212	0.069	0.021	0.090			
250155	ECX113019A	0.2979	0.2887	0.4136	0.3646	0.3990	0.3428	2.1066	6	0.351	0.051	0.160	0.045	YES
250156	ECX114019A	0.2564	0.3583	0.1491	0.2199	0.3206	0.2450	1.5493	6	0.258	0.074	0.067	0.045	YES
250157	ECX115019A	0.2433						0.2433	1	0.243	0.000	0.052	0.099	
250158	ECX116019A	0.2763	0.2837	0.2944	0.2433	0.2609	0.2566	1.6152	6	0.269	0.019	0.078	0.045	YES
250159	ECX117019A	0.2218	0.2622	0.3387	0.1973	0.3069	0.2687	1.5956	6	0.266	0.052	0.074	0.045	YES
250160	ECX118019A	0.3186	0.3909	0.3503	0.4022	0.3633	0.3925	2.2178	6	0.370	0.032	0.178	0.045	YES
250161	ECX119019A	0.2006	0.2804	0.1878	0.2035	0.2931	0.2035	1.3689	6	0.228	0.046	0.037	0.045	
250162	ECX120019A	0.2036	0.2222	0.2902	0.2629	0.2306	0.3022	1.5117	6	0.252	0.039	0.061	0.045	YES
250231	ECT111019A	0.3423	0.2139	0.3055	0.2599	0.3041	0.2608	1.6865	6	0.281	0.045	0.090	0.045	YES
250232	ECT111029A	0.3185	0.4081	0.3967	0.2515	0.2727	0.2666	1.9141	6	0.319	0.068	0.128	0.045	YES
250233	ECT112019A	0.3194	0.3335	0.3530	0.3220	0.3008	0.3307	1.9594	6	0.327	0.017	0.175	0.045	YES
250234	ECT113019A	0.1735	0.2883	0.2551	0.2080	0.2646	0.3338	1.5233	6	0.254	0.057	0.062	0.045	YES
250235	ECT114019A	0.3373	0.3482	0.2575	0.2613	0.3076	0.3835	1.8954	6	0.316	0.050	0.124	0.045	YES
250236	ECT115019A	0.3062	0.1832	0.2784	0.1817	0.2230	0.2913	1.4638	6	0.244	0.055	0.053	0.045	YES
250237	ECT116019A	0.3158	0.3675	0.4035	0.1149	0.2666	0.1370	1.6053	6	0.268	0.119	0.076	0.045	YES
250238	ECT117019A	0.4856	0.3718	0.5720	0.3282	0.4491	0.4467	2.6534	6	0.442	0.086	0.251	0.045	YES
250239	ECT118019A	0.3586	0.2935	0.2093	0.2968	0.1373	1.2955	5	0.259	0.086	0.068	0.048	YES	
250240	ECT119019A	0.1312	0.1459	0.1927	0.2581	0.1546	0.8825	5	0.177	0.051	-0.015	0.048		
250241	ECT120019A	0.2273	0.1847	0.2216	0.2712	0.2726	0.1907	1.3681	6	0.228	0.048	0.037	0.045	
251336	ECT126019A	0.3420	0.4314	0.4198	0.4126	0.3768	0.3718	2.3544	6	0.492	0.034	0.201	0.045	YES
251336D	ECT126019AD	0.3162	0.3633	0.3935	0.3852	0.2948	0.3093	2.0623	6	0.344	0.042	0.152	0.045	YES
Control-1R		0.2002	0.1644	0.1641	0.2591	0.2112	0.1795	1.1785	6	0.196	0.056	0.005	0.045	
Control-2R		0.1896	0.1902	0.1749	0.1697	0.2254	0.2089	1.1587	6	0.193	0.021	0.002	0.045	
Control-3R		0.1997	0.1595	0.1896	0.2174	0.1518	0.1909	1.1089	6	0.185	0.025	-0.007	0.045	
Combined								3.4461	18	0.191	0.027	0.000	0.041	
Pooled STD		0.055												

Table 7e. SOYBEAN Total Plant Weights (g)

Sample ID	Chen ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Central	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	0.5312	0.9130	0.9486	1.1729	1.2535	1.2150	6.0342	6	1.006	0.272	-1.159	0.271	YES
250022D	ECR111019AD	0.5576	0.8647	0.6094	0.8086	0.9382	1.1605	4.9390	6	0.823	0.221	-1.341	0.271	YES
250023	ECR111029A	0.6821	1.0748	0.9548	1.1228	1.7661		5.6006	5	1.120	0.400	-1.043	0.291	YES
250024	ECR112019A	0.8372	0.8242	0.7480	0.7922	1.0861	0.8139	5.1016	6	0.850	0.120	-1.374	0.271	YES
250025	ECR113019A	0.6821		1.7164	0.8420	0.6789	0.7642	4.6836	5	0.947	0.441	-1.228	0.291	YES
250026	ECR114019A	1.0425	0.8734	1.0699	0.7725	1.7802	1.5211	7.0596	6	1.177	0.392	-0.988	0.271	YES
250027	ECR115019A	0.5869	0.9408	0.7070	0.7820	0.7729	1.9279	5.6875	6	0.948	0.492	-1.217	0.271	YES
250028	ECR116019A	0.8874	1.4775	1.1237	1.2884	1.2332	1.2908	7.3010	6	1.217	0.198	-0.948	0.271	YES
250029	ECR117019A	0.8774	1.1517	1.2265	1.1855	0.9295	1.5200	6.8906	6	1.148	0.241	-1.016	0.271	YES
250030	ECR118019A	1.0195	1.2036	1.2840	1.1655	1.1843	0.8017	6.7586	6	1.126	0.140	-1.038	0.271	YES
250031	ECR119019A	1.2682	1.5047	1.8141	1.6026	1.1882	1.1537	8.5315	6	1.422	0.262	-0.743	0.271	YES
250032	ECR120019A	0.9962	1.3383	1.5571	1.4024	1.1707	1.8364	8.3031	6	1.384	0.295	-0.781	0.271	YES
250152	ECX111019A	1.6367	1.3295	1.3745	1.6642	1.4139	1.2436	8.6624	6	1.444	0.170	-0.721	0.271	YES
250152D	ECX111019AD	0.8415	0.8999	1.2818	1.3779	0.9852	1.0454	6.4317	6	1.072	0.214	-1.093	0.271	YES
250153	ECX111029A	0.8929	1.3894	1.1079	1.1253	1.5377	1.5538	7.6670	6	1.268	0.267	-0.897	0.271	YES
250154	ECX112019A	0.9085	1.3509	0.7523	0.9180	1.0512	1.1240	6.1049	6	1.017	0.208	-1.147	0.271	YES
250155	ECX113019A	0.7416	1.3943	1.2058	1.3449	1.3306	1.2262	7.2434	6	1.207	0.239	-0.957	0.271	YES
250156	ECX114019A	1.0551	1.4902	1.6150	1.2477	1.6306	1.4504	8.4890	6	1.417	0.224	-0.750	0.271	YES
250157	ECX115019A	0.5373	0.5110	0.6438	0.2296	0.6082	1.8290	4.3589	6	0.726	0.560	-1.438	0.271	YES
250158	ECX116019A	0.7725	1.6619	0.9923	1.0447	1.5697	1.5101	7.5512	6	1.259	0.178	-0.906	0.271	YES
250159	ECX117019A	0.9034	1.6056	1.7318	1.1881	1.3966	1.7501	8.5756	6	1.429	0.335	-0.735	0.271	YES
250160	ECX118019A	1.0165	1.2329	1.2577	1.2216	1.1918	1.9599	7.8804	6	1.112	0.328	-0.851	0.271	YES
250161	ECX119019A	0.5571	0.9860	0.8893	0.8242	2.0023	0.2469	5.5058	6	0.918	0.595	-1.247	0.271	YES
250162	ECX120019A	1.2680	1.0337	1.3810	0.8932	1.2338	1.6768	7.4865	6	1.248	0.273	-0.917	0.271	YES
250231	ECT111019A	1.0407	1.4869	1.3914	1.4213	0.9699	1.2954	7.6056	6	1.268	0.214	-0.897	0.271	YES
250232	ECT111029A	1.0523	1.5220	1.0687	0.8622	1.1618	1.2770	6.9440	6	1.157	0.225	-1.007	0.271	YES
250233	ECT112019A	0.9463	1.6986	1.3405	1.3853	1.5020	1.5152	8.3884	6	1.398	0.254	-0.767	0.271	YES
250234	ECT113019A	1.1423	1.3233	1.5075	1.4046	1.6199	1.7248	8.7224	6	1.454	0.210	-0.711	0.271	YES
250235	ECT114019A	0.2170	0.7262	0.6824	0.9152	0.8421	3.3829	5	0.877	0.273	-1.488	0.291	YES	
250236	ECT115019A	0.8256	1.2587	1.0063	1.1373	1.0091		5.2370	5	1.047	0.162	-1.117	0.291	YES
250237	ECT116019A	1.0715		0.9865	1.0396		0.7989	3.8965	4	0.974	0.122	-1.191	0.318	YES
250238	ECT117019A	0.8308	1.2309	0.9794	1.4062	1.0195	1.4817	6.9485	6	1.158	0.257	-1.007	0.271	YES
250239	ECT118019A		1.3725					1.3725	1	1.373	0.009	-0.792	0.594	YES
250240	ECT119019A							0.0000	0	0.000	0.000	-2.165	NA	YES
250241	ECT120019A	1.0503	0.5884	0.8326	1.7267		4.1980		4	1.050	0.489	-1.115	0.318	YES
251336	ECT126019A	0.8718	1.1420	0.8008	1.3057	0.9706	1.1980	6.2889	6	1.048	0.198	-1.117	0.271	YES
251336D	ECT126019AD	0.8617	0.4091	1.3940	0.7882	0.4045	0.8820	4.7395	6	0.790	0.667	-1.375	0.271	YES
Control-1R		1.5071	1.4931	1.1403	0.6111	2.5954	0.9316	8.2786	6	1.180	0.686	-0.785	0.271	
Control-2R		2.7015	2.0128	3.4735	2.4333	2.4571	2.4285	15.5067	6	2.584	0.489	0.420	0.271	
Control-3R		2.6740	2.7327		2.5874	2.8550	10.8493		4	2.712	0.112	0.548	0.318	
Combined								34.6344	16	2.765	0.709	0.000	0.197	
Pooled STD		0.33												

010026

Southwest Research Institute

December 16, 2004

TABLE 8a. CORN Total Root Weights (g)

Sample ID	Cheat ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	1.1142	1.1857	1.2762	2.9318	0.9128	1.2734	8.6941	6	1.449	0.739	0.799	0.293	YES
250022D	ECR111019AD	1.3615	3.0236	3.6426	1.3929	1.9285	0.9493	12.2984	6	2.050	1.060	1.399	0.293	YES
250023	ECR111029A	0.6719	0.5451	0.8178	0.9572	0.7876	0.6506	4.4302	6	0.738	0.146	0.088	0.293	
250024	ECR112019A	1.0321	2.1085	0.5714	0.8639	0.8230	0.9523	6.3512	6	1.059	0.538	0.408	0.293	YES
250025	ECR113019A	0.5310	0.6671	0.6752	3.7093	0.6534	0.9898	7.2258	6	1.204	1.237	0.554	0.293	YES
250026	ECR114019A	0.5726	1.3779	0.6194	0.7799	1.3878	1.2841	6.0217	6	1.004	0.382	0.353	0.293	YES
250027	ECR115019A	0.9003	1.5099	1.0264	2.2719	1.2123	1.1512	8.0720	6	1.345	0.498	0.695	0.293	YES
250028	ECR116019A	0.6257	0.6218	0.5812	0.7342	0.6264	0.7085	3.8978	6	0.650	0.059	-0.001	0.293	
250029	ECR117019A	0.7525	0.7606	0.7061	0.7128	0.5104	0.8639	4.3123	6	0.719	0.117	0.068	0.293	
250030	ECR118019A	0.6798	0.5126	0.7623	0.7317	0.6836	0.5817	3.9537	6	0.659	0.064	0.008	0.293	
250031	ECR119019A	0.9796	0.7002	0.8118	0.6434	1.2139	0.6405	4.9894	6	0.832	0.227	0.181	0.293	
250032	ECR120019A	0.5766	0.7497	0.6039	0.9689	0.5158	1.072	4.5221	6	0.754	0.237	0.103	0.293	
250152	ECX111019A	0.5550	0.2613	0.5124	0.4968	0.7242	0.6284	3.1781	6	0.530	0.156	-0.121	0.293	
250152D	ECX111019AD	0.5993	0.5112	0.5401	0.5426	0.5531	0.8953	3.6416	6	0.607	0.144	-0.044	0.293	
250153	ECX111029A	0.5512	0.6912	0.6768	0.6360	0.3700	0.9988	3.9240	6	0.654	0.201	0.004	0.293	
250154	ECX112019A	0.5689	0.5236	0.6919	0.5732	0.6975	0.6223	3.6774	6	0.613	0.071	-0.038	0.293	
250155	ECX113019A	0.4796	0.9198	0.9451	1.0399	0.8299	0.7240	4.9383	6	0.821	0.200	0.173	0.293	
250156	ECX114019A	0.6082	0.8597	0.5237	0.8663	0.7954	0.5489	4.2022	6	0.700	0.158	0.050	0.293	
250157	ECX115019A	0.3121	0.3947	0.5786	0.5384	0.6224	0.7628	3.2990	6	0.535	0.162	-0.116	0.293	
250158	ECX116019A	0.6716	0.6293	0.6917	0.8102	0.5667	0.5624	3.9319	6	0.655	0.092	0.005	0.293	
250159	ECX117019A	0.8890	0.4329	0.7159	0.6738	0.8971	0.8801	4.4888	6	0.748	0.182	0.098	0.293	
250160	ECX118019A	0.5262	0.6665	0.6835	0.6142	0.7384	0.5517	3.7805	6	0.630	0.081	-0.020	0.293	
250161	ECX119019A	0.6214	0.4969	0.7870	0.7328	0.5668	0.4220	3.6269	6	0.604	0.139	-0.046	0.293	
250162	ECX120019A	0.5872	0.7314	0.6065	0.6687	0.4655	0.4822	3.5415	6	0.590	0.104	-0.060	0.293	
250231	ECT111019A	0.6850	0.7678	0.8266	0.7342	0.7339	0.9998	4.7473	6	0.791	0.112	0.141	0.293	
250232	ECT111029A	0.5238	0.3838	0.5284	0.6372	0.9052	0.6447	3.6231	6	0.604	0.176	-0.047	0.293	
250233	ECT112019A	1.0152	0.3652	0.9361	0.6543	0.3971	0.7928	4.1607	6	0.693	0.272	0.043	0.293	
250234	ECT113019A	0.6164	0.8758	0.7991	0.9770	0.7066	0.9228	4.9577	6	0.826	0.129	0.176	0.293	
250235	ECT114019A	0.2890	0.3225	0.2120	0.1776	0.4949	0.4685	1.9645	6	0.327	0.131	-0.123	0.293	YES
250236	ECT115019A							0.0000	0	0.000	0.000	-0.650	NA	YES
250237	ECT116019A		0.4772	0.3910	0.3415	0.3137		1.5234	4	0.381	0.072	0.270	0.345	
250238	ECT117019A	0.5929	0.8000		0.6476	0.9621	0.4187	3.4213	5	0.684	0.355	0.034	0.315	
250239	ECT118019A		0.5580		0.3839			0.9419	2	0.471	0.123	-0.180	0.468	
250240	ECT119019A							0.0000	0	0.000	0.000	-0.650	NA	YES
250241	ECT120019A	0.3454		0.7316				1.0770	2	0.539	0.271	-0.112	0.468	
251336	ECT126019A	0.3976	0.5616	0.3154	0.6042	0.3232	0.6284	2.8304	6	0.472	0.143	-0.179	0.293	
251336(i)	ECT126019AD	0.6962	0.3826	0.4983	0.6600	0.4989	0.2367	2.9727	6	0.493	0.172	-0.153	0.293	
Control-1R		0.5691	0.5569	0.6764	0.5048	0.4803	0.5289	3.3164	6	0.553	0.069	-0.098	0.293	
Control-2R		0.6097	0.6720	0.7600	0.4239	0.9309	0.7456	4.1421	6	0.690	0.169	0.040	0.293	
Control-3R		0.9721	0.5680	0.6165	0.8846	0.4553	0.7538	4.2503	6	0.708	0.197	0.038	0.293	
Combined								11.7088	18	0.630	0.163	0.000	0.204	
Pooled STD		0.362												

**010027**

*Southwest Research Institute*

December 16, 2004

**TABLE 8b. ONION Total Root Weights (g)**

Sample ID	Cheat ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Statistical Difference @ 95%
250022	ECR111019A	0.0259	0.0208	0.0130	0.0235	0.0277	0.0509	0.1618	6	0.027	0.013	0.015	0.009 YES
250023D	ECR111019AD	0.0248	0.0238	0.0232	0.0108	0.0349	0.0407	0.1582	6	0.026	0.010	0.014	0.009 YES
250023	ECR111029A	0.0070		0.0160	0.0133	0.0167	0.0296	0.0826	5	0.017	0.008	0.004	0.009
250024	ECR112019A	0.0277	0.0362	0.0517	0.0308	0.0201	0.0359	0.2024	6	0.034	0.011	0.021	0.009 YES
250025	ECR113019A	0.0135	0.0123	0.0174	0.0159	0.0229	0.0188	0.1008	6	0.017	0.004	0.004	0.009
250026	ECR114019A	0.0165	0.0192	0.0351	0.0223	0.0239	0.0269	0.1439	6	0.024	0.007	0.012	0.009 YES
250027	ECR115019A	0.0233	0.0386	0.0289	0.0291	0.0293	0.0345	0.1837	6	0.031	0.005	0.018	0.009 YES
250028	ECR116019A	0.0048	0.0125	0.0040	0.0389	0.025	0.0051	0.0978	6	0.016	0.015	0.004	0.009
250029	ECR117019A	0.0101		0.0066	0.0015	0.0084	0.0016	0.0282	5	0.006	0.004	-0.007	0.009
250030	ECR118019A	0.0177	0.0168	0.0167	0.0249	0.0102	0.0390	0.1253	6	0.021	0.010	0.009	0.009 YES
250031	ECR119019A	0.0047	0.0160		0.0165	0.0197		0.0569	4	0.014	0.007	0.002	0.010
250032	ECR120019A	0.0234	0.0037	0.0232	0.0101	0.0675		0.1279	5	0.026	0.025	-0.013	0.009
250152	ECX111019A	0.0049	0.0100			0.0186		0.0335	3	0.011	0.007	-0.001	0.011
250152D	ECX111019AD	0.0032	0.0355	0.0119	0.0467	0.0209	0.0058	0.1240	6	0.021	0.017	0.008	0.009
250153	ECX111029A		0.0028	0.0123		0.0051	0.0075	0.0277	4	0.007	0.004	-0.005	0.010
250154	ECX112019A	0.0157	0.0040	0.0086	0.0016	0.0074	0.0165	0.0538	6	0.009	0.006	-0.001	0.009
250155	ECX113019A	0.0038	0.0292		0.0008	0.0240	0.0224	0.0802	5	0.016	0.013	0.004	0.009
250156	ECX114019A	0.0409	0.0253	0.0390	0.0230	0.0338	0.0332	0.1952	6	0.031	0.007	0.020	0.009 YES
250157	ECX115019A		0.0084				0.0084		1	0.008	0.000	-0.004	0.019
250158	ECX116019A	0.0082	0.0249	0.0095	0.0105	0.0287	0.0322	0.1140	6	0.019	0.011	0.007	0.009
250159	ECX117019A	0.0174	0.0199	0.0182	0.0562	0.0243	0.0102	0.1462	6	0.024	0.016	0.012	0.009 YES
250160	ECX118019A	0.0114	0.0244	0.0264	0.0241	0.0226	0.0207	0.1296	6	0.022	0.005	0.009	0.009 YES
250161	ECX119019A	0.0083	0.0050	0.0092	0.0022	0.0005	0.0114	0.0366	6	0.006	0.004	-0.006	0.009
250162	ECX120019A	0.0115	0.0296	0.0401	0.0276	0.0199	0.0380	0.1667	6	0.028	0.011	0.015	0.009 YES
250231	ECT111019A		0.0026	0.0000	0.0052	0.0217	0.0179	0.0474	5	0.009	0.010	-0.003	0.009
250232	ECT111029A	0.0085		0.0119	0.0080	0.0080		0.0364	4	0.009	0.002	-0.003	0.010
250233	ECT112019A	0.0099	0.0170	0.0169	0.0115	0.0070	0.0108	0.0731	6	0.012	0.004	0.000	0.009
250234	ECT113019A	0.0176	0.0375	0.0104	0.0099	0.0155	0.0141	0.1050	6	0.018	0.010	0.005	0.009
250235	ECT114019A	0.0134	0.0045	0.0035				0.0214	3	0.007	0.003	-0.005	0.011
250236	ECT115019A		0.0165		0.0015	0.0011		0.0191	3	0.006	0.009	-0.006	0.011
250237	ECT116019A							0.0000	0	0.000	0.000	-0.012	NA YES
250238	ECT117019A	0.0068			0.0309	0.0270	0.0150	0.0797	4	0.020	0.011	0.008	0.010
250239	ECT118019A					0.0175		0.0175	1	0.018	0.000	0.005	0.019
250240	ECT119019A			0.0019				0.0019	1	0.002	0.000	-0.010	0.019
250241	ECT120019A				0.0021	0.0176		0.0197	2	0.010	0.011	-0.002	0.013
251336	ECT120019A			0.0033				0.0033	1	0.003	0.000	-0.009	0.019
251336D	ECT126019AD	0.0010	0.0013			0.0036	0.0077		3	0.005	0.001	-0.010	0.011
Control-1R		0.0128	0.0108	0.0070	0.0080	0.0093	0.0154	0.0633	6	0.011	0.001	-0.002	0.009
Control-2R		0.0082	0.0129	0.0396		0.0129		0.0736	4	0.018	0.014	0.006	0.010
Control-3R		0.0071	0.0059	0.0198	0.0037		0.0112	0.0477	5	0.010	0.006	-0.003	0.009
Combined								0.1846	15	0.012	0.009	0.000	0.006
Pooled STD		0.0102											

010028

*Southwest Research Institute*

December 16, 2004

TABLE 8c. CARROT Total Root Weights (g)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Significance
250022	ECR111019A	0.0711	0.0476	0.0747	0.1280	0.0834	0.0488	0.4536	6	0.076	0.029	0.065	0.013	YES
250022D	ECR111019AD	0.0570	0.0663	0.0457	0.0639	0.0632	0.0509	0.3470	6	0.058	0.008	0.047	0.013	YES
250023	ECR111029A	0.0233	0.0351	0.0312	0.0249	0.0103	0.0454	0.1702	6	0.028	0.012	0.018	0.013	YES
250024	ECR112019A	0.0291	0.0178	0.0709	0.0386	0.0391	0.0473	0.2428	6	0.040	0.018	0.030	0.013	YES
250025	ECR113019A	0.0453	0.0410	0.0384	0.0654	0.0525	0.0532	0.2988	6	0.049	0.010	0.039	0.013	YES
250026	ECR114019A	0.0467	0.0652	0.0474	0.0689	0.0395	0.0637	0.3314	6	0.055	0.012	0.045	0.013	YES
250027	ECR115019A	0.0300	0.0581	0.0476	0.0871	0.0519	0.0596	0.3343	6	0.056	0.019	0.045	0.013	YES
250028	ECR116019A	0.0399	0.0464	0.0417	0.0195	0.0321	0.0352	0.2148	6	0.036	0.009	0.025	0.013	YES
250029	ECR117019A	0.0066	0.0022	0.0631	0.0183	0.0182	0.0139	0.1223	6	0.020	0.022	0.010	0.013	
250030	ECR118019A	0.0441	0.0468	0.0439	0.0545	0.0532	0.0564	0.2989	6	0.050	0.006	0.039	0.013	YES
250031	ECR119019A	0.0302	0.0623	0.0476	0.0289	0.0724	0.0747	0.3161	6	0.053	0.020	0.042	0.013	YES
250032	ECR120019A							0.0000	0	0.000	0.000	-0.010	NA	YES
250152	ECX111019A	0.0243	0.0010	0.0117	0.0263	0.0020	0.0653	0.1663	5	0.013	0.012	0.003	0.014	
250152D	ECX111019AD	0.0333	0.0440	0.0158	0.0113	0.0107	0.0630	0.1783	6	0.030	0.021	0.019	0.013	YES
250153	ECX111029A	0.0227	0.0273	0.0773	0.0600	0.0067	0.0627	0.2563	6	0.043	0.028	0.032	0.013	YES
250154	ECX112019A	0.0009	0.0071	0.0011	0.0299	0.0390	0	0.010	4	0.010	0.014	-0.001	0.015	
250155	ECX113019A	0.0251	0.0427	0.0484	0.0530	0.0302	0.0389	0.2383	6	0.040	0.011	0.029	0.013	YES
250156	ECX114019A	0.0320	0.0357	0.0479	0.0312	0.0591	0.0285	0.2344	6	0.039	0.012	0.029	0.013	YES
250157	ECX115019A	0.0010	0.0026	0.0092	0.0016	0.0144	0	0.004	4	0.004	0.004	-0.007	0.015	
250158	ECX116019A	0.0435	0.0402	0.0621	0.0240	0.0228	0.0807	0.2733	6	0.046	0.022	0.035	0.013	YES
250159	ECX117019A	0.0122	0.0530	0.0332	0.0125	0.0395	0.0460	0.1964	6	0.033	0.017	0.022	0.013	YES
250160	ECX118019A	0.0366	0.0139	0.0347	0.0395	0.0435	0.0171	0.1853	6	0.031	0.012	0.020	0.013	YES
250161	ECX119019A	0.0081	0.0031	0.0079	0.0168	0.0093	0.0080	0.0532	6	0.009	0.004	-0.002	0.013	
250162	ECX120019A	0.0484	0.0637	0.0316	0.0600	0.0277	0.0294	0.2608	6	0.043	0.016	0.033	0.013	YES
250231	ECT111019A	0.0279	0.0431	0.0589	0.0452	0.0714	0.2465	0.6806	5	0.049	0.017	0.039	0.014	YES
250232	ECT111029A	0.0374	0.0415	0.0289	0.0152	0.0077	0.1307	0.4284	5	0.026	0.014	0.016	0.014	YES
250233	ECT112019A	0.0372	0.0288	0.0207	0.0165	0.0160	0.0086	0.1278	6	0.021	0.010	0.011	0.013	
250234	ECT113019A	0.0313	0.0269	0.0516	0.0235	0.0247	0.0220	0.1800	6	0.030	0.017	0.020	0.013	YES
250235	ECT114019A	0.0023	0	0.0250	0	0.0273	0	0.014	2	0.014	0.016	0.003	0.011	
250236	ECT115019A	0.0064	0	0	0	0.0064	0	0.006	1	0.006	0.000	-0.004	0.028	
250237	ECT116019A	0	0	0	0	0.0000	0	0.000	0	0.000	0.000	-0.010	NA	YES
250238	ECT117019A	0	0	0	0	0.0004	1	0.000	1	0.000	0.000	-0.010	0.028	
250239	ECT118019A	0	0	0	0	0.0000	0	0.000	0	0.000	0.000	-0.010	NA	YES
250240	ECT119019A	0	0	0	0	0.0000	0	0.000	0	0.000	0.000	-0.010	NA	YES
250241	ECT120019A	0	0	0	0	0.0025	0	0.0025	1	0.003	0.000	-0.008	0.028	
251336	ECT126019A	0	0	0.0339	0	0.0339	1	0.034	1	0.000	0.023	0.028		
251336D	ECT126019AD	0	0	0.0005	0	0.0131	0	0.0136	2	0.007	0.009	-0.004	0.021	
Control-1R		0.0300	0.0060	0.0073	0.0039	0.0039	0.0511	0.0511	5	0.010	0.011	0.000	0.014	
Control-2R		0.0167	0.0213	0.0169	0.0087	0.0127	0.0085	0.0848	6	0.014	0.005	0.004	0.013	
Control-3R		0.0077	0	0.0007	0	0.0018	0	0.0102	3	0.003	0.004	-0.007	0.017	
Combined		0	0	0	0	0	0.1461	0.1461	14	0.010	0.008	0.000	0.010	
Pooled STD		0.015												

010029

*Southwest Research Institute*

December 16, 2004

TABLE 8d. RADISH Total Root Weights (g)

Sample ID	Chen ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Central	Least Significant Difference @ 95% Statistical Difference
250022	ECR111019A	0.0876	0.2513	0.2795	0.2275	0.3184	0.2890	1.4533	6	0.242	0.082	0.152	0.054 YES
250022D	ECR111019AD	0.130	0.2519	0.2464	0.1632	0.1036	0.1898	1.2679	6	0.211	0.074	0.121	0.054 YES
250023	ECR111029A	0.1130	0.1126	0.1596	0.1302	0.1074	0.1355	0.7583	6	0.126	0.020	0.036	0.054
250024	ECR112019A	0.0999	0.1237	0.1314	0.1019	0.1230	0.0841	0.6640	6	0.111	0.018	0.020	0.054
250025	ECR113019A	0.0870	0.0925	0.1019	0.0801	0.0608	0.1293	0.5516	6	0.092	0.023	0.001	0.054
250026	ECR114019A	0.0757	0.0907	0.0532	0.0812	0.0538	0.0798	0.4344	6	0.072	0.015	-0.018	0.054
250027	ECR115019A	0.1193	0.9880	0.1415	0.1257	0.1029	0.1050	1.5824	6	0.264	0.155	0.173	0.054 YES
250028	ECR116019A	0.0629	0.0959	0.1068	0.0837	0.0750	0.0935	0.5178	6	0.086	0.016	-0.004	0.054
250029	ECR117019A	0.0854	0.0819	0.1300	0.0738	0.0500	0.0827	0.5038	6	0.084	0.026	-0.007	0.054
250030	ECR118019A	0.1035	0.1017	0.0499	0.0504	0.0726	0.0637	0.4468	6	0.074	0.024	-0.016	0.054
250031	ECR119019A	0.1545	0.1395	0.1625	0.1460	0.1245	0.1312	0.8578	6	0.143	0.014	0.052	0.054
250032	ECR120019A	0.0918	0.1292	0.1014	0.0580	0.1495	0.0558	0.5857	6	0.098	0.018	0.007	0.054
250152	ECX111019A	0.0510	0.1110		0.0692	0.1204	0.1269	0.4785	5	0.096	0.014	0.005	0.058
250152D	ECX111019AD	0.0526	0.1219	0.0670	0.0833	0.0702	0.0668	0.4618	6	0.077	0.024	-0.014	0.054
250153	ECX111029A	0.1273	0.0809	0.0829	0.0626	0.0741	0.0808	0.5084	6	0.083	0.022	-0.006	0.054
250154	ECX112019A	0.0628			0.0000	0.0150	0.0728		3	0.026	0.031	-0.065	0.071
250155	ECX113019A	0.0609	0.0633	0.0964	0.1126	0.0952	0.1272	0.5556	6	0.093	0.026	0.002	0.054
250156	ECX114019A	0.0567	0.0699	0.0442	0.0665	0.0949	0.0755	0.4077	6	0.068	0.017	-0.023	0.054
250157	ECX115019A		0.0647				0.0647		1	0.065	0.000	-0.026	0.118
250158	ECX116019A	0.0470	0.0836	0.0389	0.0587	0.0404	0.0444	0.3130	6	0.052	0.017	-0.018	0.054
250159	ECX117019A	0.0724	0.0918	0.1132	0.0995	0.1045	0.1109	0.5923	6	0.099	0.015	0.008	0.054
250160	ECX118019A	0.0793	0.0608	0.1708	0.1500	0.1333	0.1256	0.7198	6	0.120	0.042	0.029	0.054
250161	ECX119019A	0.0363	0.0338	0.0309	0.0313	0.0593	0.0338	0.2252	6	0.038	0.011	-0.053	0.054
250162	ECX120019A	0.0340	0.0441	0.0513	0.0311	0.0299	0.0331	0.2335	6	0.047	0.009	-0.053	0.054
250231	ECT111019A	0.1006	0.0942	0.1308	0.0684	0.0700	0.0796	0.5426	6	0.090	0.023	0.000	0.054
250232	ECT111029A	0.1062	0.1176	0.0987	0.0516	0.2158	0.2544	0.8443	6	0.141	0.077	0.080	0.054
250233	ECT112019A	0.0606	0.0491	0.0504	0.0501	0.0766	0.0610	0.3478	6	0.058	0.011	-0.033	0.054
250234	ECT113019A	0.0504	0.0544	0.1095	0.0778	0.1175	0.1028	0.5124	6	0.095	0.029	-0.005	0.054
250235	ECT114019A	0.1223	0.1206	0.0802	0.0755	0.0892	0.1212	0.6090	6	0.102	0.022	0.011	0.054
250236	ECT115019A	0.0622	0.0424	0.0453	0.0306	0.0304	0.0437	0.2546	6	0.042	0.017	-0.048	0.054
250237	ECT116019A	0.0614	0.0736	0.0771	0.0218	0.0477	0.0196	0.3012	6	0.050	0.025	-0.040	0.054
250238	ECT117019A	0.1311	0.1245	0.1454	0.0655	0.0933	0.1005	0.6603	6	0.110	0.029	0.019	0.054
250239	ECT118019A	0.0802	0.0662		0.0395	0.0577		0.2446	4	0.051	0.017	-0.030	0.061
250240	ECT119019A	0.0148		0.0462	0.0395	0.0376	0.0258	0.1639	5	0.033	0.012	-0.058	0.058 YES
250241	ECT120019A	0.0608	0.0769	0.0386	0.0872	0.0666	0.0248	0.3549	6	0.059	0.024	-0.031	0.054
251336	ECT126019A	0.1211	0.1380	0.1048	0.1211	0.1631	0.0472	0.6953	6	0.116	0.039	0.025	0.054
251336D	ECT126019AD	0.0698	0.1046	0.0817	0.0776	0.0736	0.0684	0.4757	6	0.079	0.013	-0.011	0.054
Control-1R		0.0927	0.0664	0.0894	0.0964	0.1265	0.1078	0.5792	6	0.097	0.020	0.006	0.054
Control-2R		0.0745	0.1966	0.1049	0.0681	0.0711	0.1120	0.6272	6	0.105	0.049	0.014	0.054
Control-3R		0.1199	0.0548	0.0612	0.0783	0.0509	0.0589	0.4240	6	0.071	0.026	-0.020	0.054
Combined								1.6304	18	0.091	0.035	0.000	0.037
Pooled STD		0.0661											

010030

Southwest Research Institute

December 16, 2004

TABLE 8e. SOYBEAN Total Root Weights (g)

Sample ID	Clean ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Central	1-sided Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	0.5312	0.9130	0.9486	1.1729	1.2535	1.2180	6.0342	6	1.006	0.272	0.127	0.114	YES
250022D	ECR111019AD	0.3033	0.3944	0.4340	0.3378	0.4142	0.5921	2.4758	6	0.413	0.101	-0.466	0.114	YES
250023	ECR111029A	0.3660	0.5148	0.5165	0.4595	0.7467		2.6035	5	0.521	0.140	-0.358	0.123	YES
250024	ECR112019A	0.1719	0.2186	0.1932	0.3248	0.3942	0.2225	1.5252	6	0.234	0.086	-0.024	0.114	YES
250025	ECR113019A	0.3760		0.6115	0.3566	0.3189	0.5294	2.2124	5	0.442	0.133	-0.436	0.123	YES
250026	ECR114019A	0.4310	0.3420	0.4251	0.6732	0.5514	0.6180	3.0407	6	0.507	0.128	-0.372	0.114	YES
250027	ECR115019A	0.2036	0.4451	0.3656	0.3843	0.6731	0.8170	2.8887	6	0.481	0.224	-0.397	0.114	YES
250028	ECR116019A	0.5607	0.5293	0.3760	0.4761	0.4259	0.3469	2.7149	6	0.452	0.085	-0.426	0.114	YES
250029	ECR117019A	0.2593	0.3443	0.3225	0.3417	0.4241	0.4336	2.1255	6	0.354	0.066	-0.524	0.114	YES
250030	ECR118019A	0.4435	0.3739	0.7295	0.4121	0.5100	0.4269	2.8959	6	0.483	0.129	-0.396	0.114	YES
250031	ECR119019A	0.6109	0.4068	0.5296	0.5358	0.3384	0.3147	2.7362	6	0.456	0.120	-0.422	0.114	YES
250032	ECR120019A	0.4399	0.3487	0.4397	0.4644	0.3212	0.5156	2.5295	6	0.422	0.073	-0.457	0.114	YES
250152	ECX111019A	0.4531	0.4223	0.4662	0.4988	0.4218	0.3688	2.6310	6	0.439	0.045	-0.440	0.114	YES
250152D	ECX111019AD	0.3501	0.2939	0.3356	0.4143	0.3362	0.2883	2.0184	6	0.336	0.046	-0.542	0.114	YES
250153	ECX111029A	0.3110	0.3619	0.3089	0.3449	0.4364	0.4094	2.1725	6	0.362	0.052	-0.516	0.114	YES
250154	ECX112019A	0.2562	0.3052	0.2284	0.2329	0.3050	0.3158	1.6435	6	0.274	0.039	-0.605	0.114	YES
250155	ECX113019A	0.2472	0.3147	0.2908	0.3216	0.3485	0.3916	1.9144	6	0.319	0.049	-0.359	0.114	YES
250156	ECX114019A	0.4633	0.4697	0.3930	0.3664	0.4477	0.4035	2.5436	6	0.424	0.042	-0.455	0.114	YES
250157	ECX115019A	0.2322	0.3003	0.3330	0.0824	0.3608	0.5630	1.8717	6	0.312	0.158	-0.567	0.114	YES
250158	ECX116019A	0.2303	0.4392	0.3616	0.3014	0.4064	0.4522	2.1911	6	0.365	0.086	-0.514	0.114	YES
250159	ECX117019A	0.2245	0.3534	0.4765	0.3422	0.4713	0.6612	2.5291	6	0.422	0.150	-0.457	0.114	YES
250160	ECX118019A	0.4922	0.4418	0.4077	0.4041	0.4041	0.5168	2.6667	6	0.444	0.049	-0.434	0.114	YES
250161	ECX119019A	0.1779	0.2943	0.2529	0.2449	0.5153	0.1308	1.6161	6	0.269	0.134	-0.609	0.114	YES
250162	ECX120019A	0.4082	0.3380	0.3694	0.2937	0.4372	0.5599	2.4064	6	0.401	0.093	-0.477	0.114	YES
250231	ECT111019A	0.4221	0.4305	0.4733	0.4201	0.3609	0.3816	2.4885	6	0.415	0.039	-0.464	0.114	YES
250232	ECT111029A	0.3017	0.4804	0.2949	0.1985	0.2807	0.3748	1.9310	6	0.322	0.096	-0.557	0.114	YES
250233	ECT112019A	0.3317	0.4325	0.3613	0.5143	0.5313	0.5283	2.6994	6	0.450	0.088	-0.429	0.114	YES
250234	ECT113019A	0.3620	0.3729	0.4967	0.3817	0.4883	0.6082	2.7098	6	0.452	0.097	-0.427	0.114	YES
250235	ECT114019A	0.0984	0.0910	0.1304	0.0913	0.1859	0.5970	5	0.119	0.047	-0.759	0.123	YES	
250236	ECT115019A	0.1980	0.4790	0.3680	0.4082	0.3987		1.8519	5	0.370	0.105	-0.508	0.123	YES
250237	ECT116019A	0.3352		0.3744	0.3268		0.3084	1.3448	4	0.336	0.028	-0.542	0.134	YES
250238	ECT117019A	0.2351	0.3655	0.2729	0.4149	0.2694	0.3644	1.9222	6	0.320	0.071	-0.558	0.114	YES
250239	ECT118019A		0.4652				0.4652		1	0.465	0.000	-0.413	0.251	YES
250240	ECT119019A						0.0000	0	0	0.000	0.000	-0.879	NA	YES
250241	ECT120019A	0.2525	0.2309	0.2777	0.4388			1.1999	4	0.300	0.093	-0.579	0.134	YES
251336	ECT126019A	0.3953	0.2898	0.3840	0.3336	0.3628	0.2793	2.0448	6	0.341	0.048	-0.538	0.114	YES
251336D	ECT126019AD	0.2260	0.2037	0.4312	0.2327	0.1415	0.2758	1.5109	6	0.252	0.059	-0.627	0.114	YES
Control-1R		0.6043	1.0614	0.7138	0.2405	1.3908	0.3458	4.3566	6	0.726	0.436	-0.152	0.114	
Control-2R		0.8311	0.6571	1.5250	0.8149	0.8620	1.1399	5.8300	6	0.972	0.313	0.093	0.114	
Control-3R		0.9012	0.8865		0.9081	1.1737		3.8695	4	0.967	0.198	0.089	0.134	
Combined							14.0561	16	0.879	0.339	0.000	0.089		
Pooled STD		0.139												

010031

Southwest Research Institute

December 16, 2004

TABLE 9a. CORN Average Plant Heights (cm)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Central	Least Significant Difference @ 95%	Statistical Difference	
250032	ECR111019A	37.5	41.8	43.2	46.0	42.3	52.5	263.3	6	44.9	3.0	-0.9	4.9		
250030	ECR111019AD	47.8	53.8	49.3	35.1	47.3	44.2	277.5	6	46.3	6.3	1.5	4.9		
250023	ECR111029A	42.5	38.5	45.5	48.0	41.8	50.8	267.1	6	44.5	4.5	-0.3	4.9		
250024	ECR112019A	43.3	43.3	44.7	41.3	38.3	42.7	253.7	6	42.3	2.2	-2.5	4.9		
250025	ECR113019A	35.5	48.2	48.8	46.5	45.7	47.0	271.6	6	45.3	4.9	0.5	4.9		
250026	ECR114019A	44.5	52.0	47.8	51.0	45.7	49.8	290.8	6	48.5	3.0	3.7	4.9		
250027	ECR115019A	43.0	44.3	52.5	43.8	47.7	48.7	280.0	6	46.7	3.6	-1.9	4.9		
250028	ECR116019A	42.5	44.8	47.0	37.8	46.3	43.0	261.5	6	45.6	3.3	-1.2	4.9		
250029	ECR117019A	33.0	35.7	41.0	41.3	48.0	50.0	249.0	6	41.3	6.6	-3.3	4.9		
250030	ECR118019A	41.2	44.8	49.7	51.0	47.8	54.0	288.5	6	48.1	4.6	3.3	4.9		
250031	ECR119019A	33.0	28.0	47.3	49.0	47.8	51.5	256.7	6	42.8	9.7	-2.0	4.9		
250032	ECR120019A	34.8	29.7	46.8	40.3	41.0	47.5	240.1	6	40.0	6.9	-4.8	4.9		
250152	ECX111019A	23.8	38.8	41.0	38.5	38.2	34.5	214.8	6	35.8	6.2	-9.0	4.9	YES	
250152D	ECX111019AD	41.7	42.3	28.8	41.7	49.2	42.0	247.7	6	41.3	6.7	-3.5	4.9		
250153	ECX111029A	28.3	43.2	44.3	46.3	26.5	46.7	235.3	6	39.2	9.3	-5.6	4.9	YES	
250154	ECX112019A	39.2	45.3	45.0	36.5	41.5	41.3	248.5	6	41.3	6.4	-3.4	4.9		
250155	ECX113019A	44.5	29.2	47.5	51.5	40.3	51.0	264.0	6	44.0	8.4	-0.8	4.9		
250156	ECX114019A	34.8	41.2	50.3	46.2	47.5	52.5	272.4	6	45.4	6.5	0.6	4.9		
250157	ECX115019A	44.3	46.8	46.7	48.3	55.3	43.8	285.0	6	47.5	4.1	2.7	4.9		
250158	ECX116019A	33.7	43.7	34.5	46.7	47.8	46.8	253.0	6	42.2	6.4	-2.6	4.9		
250159	ECX117019A	51.8	54.7	47.7	47.5	48.3	43.3	293.3	6	48.9	3.9	4.1	4.9		
250160	ECX118019A	46.7	54.0	50.3	51.3	50.3	47.3	300.0	6	50.0	2.7	5.2	4.9	YES	
250161	ECX119019A	42.5	38.8	50.2	51.2	46.8	54.3	283.7	6	47.3	5.8	2.5	4.9		
250162	ECX120019A	47.0	49.7	64.5	49.0	54.5	314.3	514.3	6	52.4	6.4	7.6	4.9	YES	
250234	ECT111019A	42.7	42.3	49.7	43.7	30.0	41.3	249.7	6	41.6	6.4	-3.2	4.9		
250232	ECT111029A	38.5	44.3	42.3	38.3	44.0	35.8	243.3	6	40.5	3.5	-4.1	4.9		
250233	ECT112019A	43.3	51.5	62.0	56.8	53.0	59.2	325.8	6	54.4	6.6	9.5	4.9	YES	
250234	ECT113019A	60.8	52.7	52.8	57.8	52.3	54.5	330.9	6	53.2	3.4	10.3	4.9	YES	
250235	ECT114019A	53.5	50.0	47.0	39.0	35.7	33.3	258.5	6	43.3	8.2	-1.7	4.9		
250236	ECT115019A								0	0	0.0	0.0	-44.8	NA	YES
250237	ECT116019A		35.8	43.4	28.9	39.0		147.1	4	36.8	6.1	-8.0	5.7	YES	
250238	ECT117019A	39.7	48.7		54.0	53.0	44.0	239.3	5	47.9	6.1	3.1	5.2		
250239	ECT118019A		52.0			49.0		101.0	2	50.5	2.1	5.7	7.8		
250240	ECT119019A								0	0	0.0	0.0	-44.8	NA	YES
250241	ECT120019A	40.0		45.0				85.0	2	42.5	3.5	-2.3	7.8		
251336	ECT126019A	41.9	53.5	52.7	46.9	52.2	46.0	293.0	6	48.8	4.6	4.6	4.9		
251336D	ECT126019AD	42.6	62.0	50.6	43.2	40.4	46.0	284.7	6	47.5	7.9	2.6	4.9		
Control-1R		40.5	45.8	45.0	40.3	33.5	39.3	244.3	6	40.7	4.4	-4.1	4.9		
Control-2R		63.0	51.3	49.7	42.0	48.3	45.0	299.3	6	49.9	7.2	5.1	4.9		
Control-3R		47.0	50.0	46.3	46.0	23.5	50.0	262.8	6	43.8	10.1	-1.0	4.9		
Combined								806.5	18	44.8	8.16	0.0	3.4		
Pooled STD		6.01													

010032

Southwest Research Institute

December 16, 2004

TABLE 9b. ONION Average Plant Heights (cm)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95% Statistical Difference	
250022	ECR111019A	11.6	12.0	12.0	12.5	10.3	13.2	71.6	6	11.9	1.0	3.1	3.0	YES
250022D	ECR111019AD	11.8	12.0	13.0	15.0	14.1	16.1	82.0	6	13.7	1.7	4.9	3.0	YES
250023	ECR111029A	7.8		15.0	11.3	14.3	12.5	69.8	5	12.2	2.9	3.4	3.2	YES
250024	ECR112019A	15.5	13.0	14.0	16.0	17.0	16.3	91.8	6	15.3	1.5	6.5	3.0	YES
250025	ECR113019A	12.5	10.5	10.5	13.6	15.3	11.0	73.5	6	12.2	2.0	3.4	3.0	YES
250026	ECR114019A	11.2	14.3	13.4	14.6	15.7	14.3	83.3	6	13.9	1.5	5.1	3.0	YES
250027	ECR115019A	13.5	12.9	12.2	15.4	17.3	12.8	84.0	6	14.0	1.9	5.2	3.0	YES
250028	ECR116019A	6.5	20.5	12.0	17.1	16.3	7.5	80.0	6	14.3	5.6	4.3	3.0	YES
250029	ECR117019A	14.5		9.7	16.5	11.0	7.8	59.4	5	11.9	3.6	3.1	3.2	
250030	ECR118019A	10.4	13.8	10.7	18.0	11.0	13.9	77.8	6	13.9	2.9	4.2	3.0	YES
250031	ECR119019A	18.0	19.3	0.0	20.0	20.5	0.0	77.8	6	13.9	10.1	4.2	3.0	YES
250032	ECR120019A	16.8	12.0	16.0	15.3	16.1		76.1	5	15.2	1.9	6.4	3.2	YES
250152	ECX111019A	16.8	16.5	6.0	1.0	16.3	0.0	56.8	6	9.5	8.1	0.7	3.0	
250152D	ECX111019AD	14.5	15.9	13.5	16.9	16.8	15.0	92.6	6	15.4	1.3	6.6	3.0	YES
250153	ECX111029A	0.0	7.3	9.5		9.3	9.2	35.2	5	7.0	4.0	1.8	3.2	
250154	ECX112019A	8.5	7.5		11.5	6.5	12.8	46.8	5	9.4	2.7	0.5	3.2	
250155	ECX113019A	12.0	15.7	0.0	11.0	17.0	17.5	73.2	6	12.2	6.5	3.4	3.0	YES
250156	ECX114019A	14.5	12.9	13.5	14.0	15.0	14.8	84.8	6	14.1	0.8	3.5	3.0	YES
250157	ECX115019A			9.3				9.3	1	9.3	0.0	0.4	6.5	
250158	ECX116019A	11.8	13.3	13.0	15.8	15.8	16.4	86.0	6	14.3	1.9	5.5	3.0	YES
250159	ECX117019A	18.3	13.2	15.5	17.5	15.1	14.3	93.8	6	15.6	1.9	6.8	3.0	YES
250160	ECX118019A	13.5	16.0	13.0	17.6	16.0	15.1	91.2	6	15.2	1.7	6.4	3.0	YES
250161	ECX119019A	11.3	11.7	13.6	10.1	14.5	14.9	76.1	6	12.7	1.9	3.9	3.0	YES
250162	ECX120019A	14.9	16.3	14.6	16.9	14.3	16.6	93.7	6	15.6	1.1	6.8	3.0	YES
250231	ECT111019A		6.0		17.0	16.0	18.3	57.3	4	14.3	5.6	3.5	3.5	
250232	ECT111029A	17.0		15.3	12.5	16.5		61.3	4	15.3	2.0	6.5	3.5	YES
250233	ECT112019A	16.8	13.0	14.9	18.5	15.5	14.6	93.2	6	15.5	1.9	6.7	3.0	YES
250234	FCT113019A	12.6	13.8	8.5	16.3	18.0	15.8	85.0	6	14.2	3.4	5.4	3.0	YES
250235	ECT114019A		13.3	14.0	9.5			36.8	3	12.3	2.4	3.4	4.0	
250236	ECT115019A	15.4			15.6	3.5		34.5	3	11.5	6.9	2.7	4.0	
250237	ECT116019A							0.0	0	0.0	0.0	-8.8	N/A	YES
250238	ECT117019A	18.3			19.2	17.0	19.5	73.9	4	18.5	1.1	9.7	3.5	YES
250239	ECT118019A						10.5		1	10.5	0.0	1.7	6.5	
250240	ECT119019A				8.7			8.7	1	8.7	0.0	-0.1	6.5	
250241	ECT120019A					9.5	13.0	22.5	2	11.3	2.5	2.4	4.7	
251336	ECT126019A				11.7			11.7	1	11.7	0.0	2.9	6.5	
251336D	ECT126019AD	11.0			12.5			14.0	3	12.5	1.5	3.7	4.0	
Control-1R		9.0	8.0	7.3	6.8	7.3	10.2	48.4	6	8.1	1.3	-0.7	3.0	
Control-2R		6.0	7.8	9.6		13.0		36.4	4	9.1	3.0	0.3	3.5	
Control-3R		12.0	10.0	9.0	7.0		9.3	47.3	5	9.5	1.8	0.6	3.2	
Combined								132.0	15	8.80	1.96	0.0	2.2	
Pooled STD		3.60												

010033

Southwest Research Institute

December 16, 2004

TABLE 9c. CARROT Average Plant Heights (cm)

Sample ID	Clean ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	6.8	7.5	7.8	9.4	9.4	8.5	49.4	6	8.2	1.1	4.3	2.0	YES
250022D	ECR111019AD	8.6	7.1	6.8	9.6	8.9	8.6	49.7	6	8.3	1.1	4.4	2.0	YES
250023	ECR111029A	8.5	4.6	7.6	8.4	5.9	6.9	41.9	6	7.0	1.5	3.0	2.0	YES
250024	ECR112019A	7.4	8.6	9.0	9.5	8.0	10.3	52.8	6	8.8	1.0	4.8	2.0	YES
250025	ECR113019A	9.1	7.9	7.0	7.8	9.0	8.4	49.2	6	8.2	0.8	4.2	2.0	YES
250026	ECR114019A	7.6	7.6	7.2	8.1	7.5	6.8	44.8	6	7.5	0.5	3.5	2.0	YES
250027	ECR115019A	8.3	10.1	8.3	10.3	9.1	7.6	53.8	6	9.0	1.1	5.0	2.0	YES
250028	ECR116019A	12.5	12.2	11.0	9.1	10.6	9.6	65.0	6	10.8	1.4	6.9	2.0	YES
250029	ECR117019A	7.5	7.8	13.8	11.3	8.6	9.8	58.8	6	9.8	2.4	5.8	2.0	YES
250030	ECR118019A	9.4	10.3	10.8	11.5	9.8	9.9	61.7	6	10.3	0.8	6.4	2.0	YES
250031	ECR119019A	18.0	11.8	12.0	12.0	9.9	12.8	76.5	6	12.8	2.7	8.8	2.0	YES
250032	ECR120019A								0	0.0	0.0	-4.0	NA	YES
250152	ECX111019A	0.0	9.4	6.7	13.0	9.3	7.5	45.8	6	7.6	4.3	3.7	2.0	YES
250152D	ECX111019AD	9.5	8.8	10.3	8.1	7.3	8.5	52.6	6	8.8	1.1	4.8	2.0	YES
250153	ECX111029A	10.0	8.6	9.9	30.8	7.0	8.8	75.1	6	12.5	9.0	8.6	2.0	YES
250154	ECX112019A	8.8				8.5		5.5	3	7.6	1.8	3.6	2.6	YES
250155	ECX113019A	14.5	11.8	12.4	13.0	11.1	14.3	77.1	6	12.8	1.1	8.9	2.0	YES
250156	ECX114019A	7.7	9.1	8.8	7.8	9.3	7.6	50.3	6	8.4	0.8	4.4	2.0	YES
250157	ECX115019A		4.0	5.5		8.0	4.0	21.5	4	5.4	1.9	1.4	2.8	
250158	ECX116019A	10.4	8.9	11.8	12.5	9.8	10.4	63.8	6	10.6	1.3	6.7	2.0	YES
250159	ECX117019A	7.5	10.3	10.3	8.2	12.2	9.3	57.7	6	9.6	1.7	5.7	2.0	YES
250160	ECX118019A	9.1	10.5	10.4	8.3	10.7	9.8	58.7	6	9.8	1.0	5.8	2.0	YES
250161	ECX119019A	7.6	8.8	6.9	8.8	9.0	6.9	48.0	6	8.0	1.0	4.6	2.0	YES
250162	ECX120019A	9.4	8.3	10.8	11.1	10.3	9.3	59.1	6	9.9	1.0	5.9	2.0	YES
250231	ECT111019A	9.5	13.3	11.1		9.7	11.1	54.7	5	10.9	1.5	7.0	2.1	YES
250232	ECT111029A		8.8	8.8	11.3	10.0	11.0	49.8	5	10.0	1.2	6.0	2.1	YES
250233	ECT112019A	8.9	7.9	7.3	7.8	9.3	8.3	49.4	6	8.2	0.7	4.3	2.0	YES
250234	ECT113019A	9.1	8.8	10.5	8.8	9.3	9.9	56.3	6	9.4	0.7	5.4	2.0	YES
250235	ECT114019A		6.5			12.5		19.0	2	9.5	4.2	5.5	3.1	YES
250236	ECT115019A	7.5						7.5	1	7.5	0.0	3.5	4.3	
250237	ECT116019A							0.0	0	0.0	0.0	-4.0	NA	YES
250238	ECT117019A						5.0	5.0	1	5.0	2.0	1.0	4.3	
250239	ECT118019A						1.0	1.0	1	1.0	0.0	-3.0	4.3	
250240	ECT119019A							0.0	0	0.0	0.0	-4.0	NA	YES
250241	ECT120019A						6.5	6.5	1	6.5	0.0	2.5	4.3	
251336	ECT126019A				18.1			18.1	1	18.1	0.0	14.1	4.3	YES
251336D	ECT126019AD				5.0	15.0		20.0	2	10.0	7.1	6.0	3.1	YES
Control-1R		4.4	3.2	3.8	5.5	3.0		19.9	5	4.0	1.0	0.0	2.1	
Control-2R		3.7	4.2	3.5	4.8	4.5	5.3	26.0	6	4.3	0.7	0.4	2.0	
Control-3R		4.0		2.0			3.5	9.5	3	3.2	1.0	-0.8	2.6	
Combined								55.3	14	3.95	0.93	0.0	1.5	
Pooled STD								2.36						

010034

Southwest Research Institute

December 16, 2004

TABLE 9d. RADISH Average Plant Heights (cm)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95% Statistical Difference	
250022	ECR111019A	4.5	4.6	4.3	4.4	4.5	4.8	27.1	6	4.5	0.2	0.4	0.6	
25002D	ECR111019AD	5.5	3.9	5.0	5.1	4.3	4.8	28.6	6	4.8	0.6	0.5	0.6	
250023	ECR111029A	5.0	5.0	4.9	5.0	4.6	5.0	29.5	6	4.9	0.2	0.7	0.6	
250024	ECR112019A	4.9	4.5	4.3	5.2	4.8	4.5	28.2	6	4.7	0.1	0.5	0.6	
250025	ECR113019A	4.8	5.0	4.8	4.6	4.8	5.0	29.0	6	4.8	0.2	0.6	0.6	
250026	ECR114019A	4.5	4.5	5.1	5.0	4.9	4.6	28.6	6	4.8	0.3	0.5	0.6	
250027	ECR115019A	4.4	3.9	3.6	4.9	4.2	3.9	24.9	6	4.2	0.5	-0.1	0.6	
250028	ECR116019A	5.3	5.0	5.8	5.5	5.0	4.9	31.4	6	5.2	0.3	1.0	0.6	
250029	ECR117019A	5.7	7.0	6.5	6.8	7.0	4.6	37.6	6	6.3	1.0	2.0	0.6	
250030	ECR118019A	5.4	5.2	5.3	5.7	5.0	5.4	32.0	6	5.4	0.2	1.1	0.6	
250031	ECR119019A	8.2	5.5	6.1	6.6	5.9	5.4	37.7	6	6.4	1.0	2.0	0.6	
250032	ECR120019A	6.3	6.0	6.2	5.8	6.0	6.0	36.3	6	6.0	0.2	1.8	0.6	
250152	ECX111019A	12.5	6.5		7.5	6.6	6.8	39.9	5	8.0	2.6	3.7	0.7	
250152D	ECX111019AD	5.6	5.6	5.3	5.9	5.4	5.9	33.7	6	5.6	0.2	1.4	0.6	
250153	ECX111029A	6.3	5.9	7.1	6.1	5.4	5.5	36.3	6	6.1	0.6	1.8	0.6	
250154	ECX112019A	3.5	6.3					7.5	22.3	3	7.4	1.1	3.2	0.8
250155	ECX113019A	8.3	8.2	9.5	8.8	7.8	6.9	49.4	6	8.2	0.9	4.6	0.6	
250156	ECX114019A	4.9	4.8	4.6	5.1	4.6	4.1	28.1	6	4.7	0.3	0.4	0.6	
250157	ECX115019A							5.5	1	5.5	0.0	1.3	1.4	
250158	ECX116019A	6.7	5.3	5.6	4.9	4.8	4.9	32.2	6	5.4	0.7	1.1	0.6	
250159	ECX117019A	6.5	5.8	6.5	5.3	4.4	5.0	33.5	6	5.6	0.8	1.5	0.6	
250160	ECX118019A	8.1	7.5	7.5	10.6	7.6	8.0	49.3	6	8.2	1.2	4.0	0.6	
250161	ECX119019A	4.2	6.1	5.6	4.9	6.1	5.0	31.9	6	5.3	0.8	1.1	0.6	
250162	ECX120019A	5.9	5.1	6.5	7.0	5.0	6.1	35.6	6	5.9	0.8	1.7	0.6	
250231	ECT111019A	7.0	5.9	7.0	6.6	6.9	5.6	39.0	6	6.5	0.6	2.3	0.6	
250232	ECT111029A	7.1	5.9	6.3	6.0	6.0	6.7	38.0	6	6.3	0.5	2.1	0.6	
250233	ECT112019A	8.1	6.6	7.7	7.1	6.6	6.4	42.5	6	7.1	0.7	2.8	0.6	
250234	ECT113019A	6.4	6.2	6.5	6.3	6.1	7.1	38.6	6	6.4	0.4	2.2	0.6	
250235	ECT114019A	5.8	4.7	4.4	4.6	5.5	5.3	30.3	6	5.1	0.6	0.8	0.6	
250236	ECT115019A	5.4	5.4	4.9	5.2	4.7	5.2	30.7	6	5.1	0.3	0.9	0.6	
250237	ECT116019A	6.0	5.1	5.7	6.5	5.4	6.5	35.1	6	5.9	0.6	1.6	0.6	
250238	ECT117019A	6.1	6.5	7.6	8.8	7.9	5.7	42.6	6	7.1	1.2	2.9	0.6	
250239	ECT118019A	4.8	4.6		5.3	4.8	6.0	25.5	5	5.1	0.6	0.8	0.7	
250240	ECT119019A	6.0		6.6	7.5	7.7	6.2	34.0	5	6.8	0.8	2.6	0.7	
250241	ECT120019A	9.0	5.5	4.3	5.0	4.5	6.5	34.7	6	5.8	1.8	1.5	0.6	
251336	ECT126019A	6.1	6.2	6.0	7.5	7.0	6.6	39.5	6	6.6	0.6	2.3	0.6	
251336D	ECT126019AD	5.2	5.5	5.7	6.0	5.5	5.4	33.3	6	5.5	0.3	1.3	0.6	
Control-1R		4.3	3.9	3.9	4.4	4.3	3.9	24.7	6	4.1	0.2	-0.1	0.6	
Control-2R		4.8	4.4	4.2	3.8	4.1	4.7	26.0	6	4.3	0.4	0.1	0.6	
Control-3R		4.3	3.7	4.0	6.0	3.8	4.0	25.8	6	4.3	0.9	0.0	0.6	
Combined								76.4	18	4.25	0.54	0.0	0.4	
Pooled STD		0.73												

**010035**

Southwest Research Institute

December 16, 2004

**TABLE 9e. SOYBEAN Average Plant Heights (cm)**

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical difference
250022	ECR111019A	15.3	16.0	25.2	20.5	20.0	20.5	117.5	6	19.6	3.6	-16.0	3.3	YES
250022D	ECR111019AD	13.8	15.5	15.0	15.3	17.8	16.8	94.3	6	15.7	1.4	-19.8	3.3	YES
250023	ECR111029A	11.3	18.3	24.0	24.7	23.0		101.3	5	20.1	5.6	-15.3	3.3	YES
250024	ECR112019A	11.7	10.7	10.7	16.7	19.3	10.3	79.3	6	13.2	3.8	-22.3	3.3	YES
250025	ECR113019A	17.7		34.2	21.5	23.8	17.0	116.1	5	23.2	6.9	-12.3	3.3	YES
250026	ECR114019A	16.5	20.7	20.3	18.5	20.8	19.8	116.7	6	19.4	1.7	-16.1	3.3	YES
250027	ECR115019A	11.0	16.3	15.8	13.8	19.3	19.5	95.8	6	16.0	3.2	-19.6	3.3	YES
250028	ECR116019A	16.8	19.0	17.0	18.2	16.7	13.2	100.8	6	16.8	2.0	-18.8	3.3	YES
250029	ECR117019A	12.7	14.5	15.7	15.7	15.2	18.3	92.0	6	15.3	1.8	-20.2	3.3	YES
250030	ECR118019A	14.7	18.7	16.7	18.5	16.7	16.8	102.0	6	17.0	1.5	-18.5	3.3	YES
250031	ECR119019A	16.5	17.5	17.5	19.2	18.5	13.5	102.7	6	17.1	2.0	-18.4	3.3	YES
250032	ECR120019A	13.3	16.5	16.3	16.7	17.2	19.0	99.0	6	16.5	1.8	-19.0	3.3	YES
250152	ECX111019A	16.5	17.0	19.3	17.2	14.2	17.0	101.1	6	16.8	1.6	-18.7	3.3	YES
250152D	ECX111019AD	14.2	15.5	18.0	16.3	16.5	17.0	97.5	6	16.3	1.4	-19.3	3.3	YES
250153	ECX111029A	18.0	21.8	18.7	14.7	18.8	15.7	107.7	6	17.9	2.5	-17.6	3.3	YES
250154	ECX112019A	14.2	15.3	15.0	8.0	10.7	7.0	70.2	6	11.7	3.7	-24.8	3.3	YES
250155	ECX113019A	20.0	18.7	17.3	20.5	17.7	13.2	107.3	6	17.9	2.6	-17.7	3.3	YES
250156	ECX114019A	17.3	18.2	18.2	18.3	17.0	15.2	104.1	6	17.3	1.2	-18.2	3.3	YES
250157	ECX115019A	23.8	21.0	24.0	14.5	20.8	21.8	125.8	6	21.0	3.5	-14.6	3.3	YES
250158	ECX116019A	14.3	17.7	16.5	16.8	17.8	14.5	97.7	6	16.3	1.5	-19.4	3.3	YES
250159	ECX117019A	13.5	26.5	20.2	22.5	19.7	16.5	118.8	6	19.8	4.5	-15.7	3.3	YES
250160	ECX118019A	10.5	26.7	33.0	32.7	27.2	24.3	174.3	6	29.1	3.5	-6.5	3.3	YES
250161	ECX119019A	18.5	22.5	21.7	19.8	24.5	15.0	121.8	6	20.3	3.3	-15.2	3.3	YES
250162	ECX120019A	18.2	17.8	19.0	20.3	18.8	19.5	113.7	6	18.9	0.9	-16.6	3.3	YES
250231	ECT111019A	16.8	21.0	18.3	17.7	15.3	18.8	108.0	6	18.0	1.9	-17.5	3.3	YES
250232	ECT111029A	18.7	19.8	17.7	16.3	13.8	14.3	100.7	6	16.8	2.4	-18.8	3.3	YES
250233	ECT112019A	17.8	21.3	27.5	22.3	23.0	17.5	129.5	6	21.0	3.7	-14.0	3.3	YES
250234	ECT113019A	23.2	25.5	23.8	20.2	18.3	19.0	130.0	6	21.7	2.9	-13.9	3.3	YES
250235	ECT114019A	15.0	9.0	10.0	12.8		8.3	55.2	5	11.0	2.8	-24.5	3.3	YES
250236	ECT115019A	23.0	24.8	24.4	22.2	48.0		142.3	5	28.5	11.0	-7.1	3.6	YES
250237	ECT116019A	17.9		18.5	21.0		17.9	75.3	4	18.8	1.5	-16.7	3.9	YES
250238	ECT117019A	18.3	16.5	11.8	15.8	19.0	15.8	97.3	6	16.2	2.5	-19.1	3.3	YES
250239	ECT118019A		23.0					23.0	1	23.0	0.0	-12.5	7.2	YES
250240	ECT119019A							0.0	0	0.0	0.0	-35.5	N/A	YES
250241	ECT120019A	14.5	20.5	21.5	17.8			74.3	4	18.6	3.1	-17.0	3.9	YES
251336	ECT126019A	19.8	23.0	20.5	23.2	21.5	18.4	126.3	6	21.1	1.9	-14.5	3.3	YES
251336D	ECT126019AD	23.6	20.5	17.8	21.6	19.6	20.0	123.1	6	20.5	1.9	-15.0	3.3	YES
Control-1R		29.3	31.0	23.3	19.5	43.3	17.5	164.0	6	27.1	9.5	-8.2	3.3	
Control-2R		33.0	36.3	52.7	42.0	33.3	33.2	230.5	6	38.4	7.8	2.9	3.3	
Control-3R		48.3	55.0		33.3	37.5		174.2	4	44.5	9.9	8.0	3.9	
Combined								568.7	16	35.54	10.82	0.0	2.4	
Pooled STD		4.01												

010036

December 16, 2004

Southwest Research Institute

TABLE 10a. CORN Average Root Lengths (cm)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Last Significant Difference @ 95% Statistical
250022	ECR11019A	18.8	19.3	20.2	27.8	19.5	23.3	128.9	6	21.5	3.5	-3.4	4.0
250022D	ECR11019AD	26.7	15.8	21.3	16.0	27.8	16.2	123.8	6	20.6	5.5	-4.3	4.0 YES
250023	ECR11029A	30.0	46.0	23.0	26.3	35.8	22.3	183.4	6	30.6	9.1	5.7	4.0 YES
250024	ECR112019A	13.3	11.7	10.7	18.3	22.0	24.7	100.7	6	16.8	5.8	-8.1	4.0 YES
250025	ECR113019A	16.3	16.3	18.8	15.5	34.8	24.7	126.3	6	21.1	7.5	-3.9	4.0
250026	ECR114019A	15.5	19.5	29.2	17.8	27.8	24.7	134.4	6	22.4	5.6	-2.5	4.0
250027	ECR115019A	17.2	13.3	17.8	26.3	23.3	24.2	122.2	6	20.4	5.0	-4.6	4.0 YES
250028	LCR116019A	10.8	23.5	23.0	19.8	19.7	15.8	112.5	6	18.8	4.8	-6.2	4.0 YES
250029	ECR117019A	15.7	17.0	15.0	16.0	11.0	20.0	94.7	6	15.8	2.9	-9.1	4.0 YES
250030	ECR118019A	20.7	21.3	20.0	16.3	23.3	21.5	123.2	6	20.5	2.3	-4.4	4.0 YES
250031	ECR119019A	15.2	13.7	16.3	21.0	19.0	20.3	105.4	6	17.6	3.0	-7.3	4.0 YES
250032	LCR120019A	15.7	18.3	17.0	32.0	24.3	17.8	125.1	6	20.8	6.2	-6.1	4.0 YES
250152	ECX111019A	21.8	21.0	23.8	19.3	20.0	27.7	133.6	6	22.3	3.1	-2.7	4.0
250152D	ECX111019AD	21.0	18.0	20.3	20.8	16.7	28.0	124.8	6	20.8	3.9	-4.1	4.0 YES
250153	ECX111029A	21.0	13.2	18.0	17.5	19.7	29.2	118.5	6	19.8	5.4	-5.2	4.0 YES
250154	ECX112019A	18.8	13.3	21.5	19.2	23.2	15.7	111.6	6	18.6	3.7	-6.3	4.0 YES
250155	ECX113019A	15.3	13.0	28.0	33.5	25.3	17.2	132.3	6	22.1	8.1	-2.9	4.0
250156	ECX114019A	23.0	13.2	19.3	19.7	20.5	26.3	121.9	6	20.3	4.4	-4.6	4.0 YES
250157	ECX115019A	16.5	16.0	22.2	14.8	15.8	25.2	110.4	6	18.4	4.2	-6.5	4.0 YES
250158	ECX116019A	13.7	22.5	20.5	15.5	21.5	12.8	106.4	6	17.7	4.3	-7.2	4.0 YES
250159	ECX117019A	15.3	20.2	18.0	22.7	18.5	14.0	108.7	6	18.1	3.2	-6.8	4.0 YES
250160	ECX118019A	15.3	11.3	15.3	17.2	20.0	12.8	91.9	6	15.3	3.1	-9.6	4.0 YES
250161	ECX119019A	22.2	11.2	16.5	20.0	19.3	12.5	101.7	6	16.9	4.4	-8.0	4.0 YES
250162	ECX120019A	13.0	19.3	21.0	8.0	31.5	17.0	109.8	6	18.1	8.0	-6.6	4.0 YES
250231	ECT111019A	15.7	22.7	21.8	19.0	14.8	17.7	111.7	6	18.6	3.2	-6.3	4.0 YES
250232	ECT111029A	16.0	16.5	21.3	20.3	23.5	22.0	119.7	6	19.9	3.0	-5.0	4.0 YES
250233	ECT112019A	17.0	20.3	9.0	17.0	22.0	19.0	104.3	6	17.4	4.5	-7.5	4.0 YES
250234	ECT113019A	20.8	26.7	18.3	25.3	27.8	17.0	135.9	6	22.7	4.6	-2.3	4.0
250235	ECT114019A	13.0	9.0	18.0	15.0	15.3	14.8	85.2	6	14.2	3.6	-10.7	4.0 YES
250236	ECT115019A							0.0	0	0.0	0.0	-24.9	NA YES
250237	ECT116019A	14.6	18.4	10.3	9.3			52.6	4	13.1	4.2	-11.8	4.7 YES
250238	EHT117019A	15.0	15.2		8.0	18.5	15.0	71.7	5	14.3	3.3	-10.6	4.3 YES
250239	EHT118019A		10.0			34.0		44.0	2	22.0	17.0	-2.9	6.3
250240	EHT119019A	9.0		15.3				0.0	0	0.0	0.0	-24.9	NA YES
250241	EHT120019A							24.3	2	12.2	4.5	-12.8	6.3 YES
251336	EHT126019A	24.8	15.2	20.9	26.0	18.6	20.6	126.0	6	21.0	4.0	-3.9	4.0
251336D	EHT126019AD	22.8	13.8	23.2	20.2	21.8	24.0	125.7	6	21.0	3.7	-4.0	4.0
Control-1R		29.0	28.5	31.5	28.5	32.5	27.0	177.0	6	29.5	2.1	4.6	4.0
Control-2R		24.0	29.7	28.3	19.0	28.7	28.5	158.2	6	26.4	4.1	1.4	4.0
Control-3R		20.5	23.0	19.3	19.0	13.5	18.0	113.3	6	18.9	3.1	-6.0	4.0
Combined								448.5	18	24.9	5.5	0.0	2.8
Pooled STD								4.91					

010037

December 16, 2004

Southwest Research Institute

TABLE 10b. ONION Average Root Lengths (cm)

Sample ID	Client ID	Container						Total	Number of Measurements	Standard Deviation	Difference from Control	Least Significant Difference @ 95% Statistical Difference
		1	2	3	4	5	6					
250022	ECR111019A	16.8	13.9	22.8	15.5	15.0	16.6	100.5	6	16.7	3.1	4.4
250022D	ECR111019AD	14.0	15.0	16.3	19.0	11.9	17.4	93.6	6	15.6	2.5	3.2
250023	ECR111029A	11.0		21.5	21.0	19.3	13.5	86.3	5	17.1	4.7	4.9
250024	ECR112019A	6.5	9.7	10.6	11.0	10.0	11.0	58.8	6	9.8	1.7	-2.6
250025	ECR113019A	14.5	19.8	14.5	13.5	18.5	14.2	94.9	6	15.8	2.6	3.5
250026	ECR114019A	9.8	19.0	19.1	17.4	15.2	12.6	91.1	6	15.5	3.7	3.2
250027	ECR115019A	13.3	13.0	17.6	12.9	14.8	13.4	84.9	6	14.1	1.8	1.8
250028	ECR116019A	6.7	18.5	9.0	9.4	16.7	2.8	63.0	6	10.5	6.0	-1.9
250029	ECR117019A	2.3		2.8	3.5	4.2	3.5	16.3	5	3.3	0.7	-9.1
250030	ECR118019A	6.5	9.7	6.8	13.6	7.7	15.3	59.5	6	9.9	3.7	-2.4
250031	ECR119019A	6.5	19.8	0.0	12.0	21.8		60.0	5	12.0	9.1	-0.4
250032	ECR120019A	10.5	7.0	13.3	8.0	24.1		62.9	5	12.6	6.9	0.2
250152	ECX111019A	7.3	5.7	1.0		8.7		22.6	4	5.6	3.1	-6.7
250152D	ECX111019AD	8.0	17.5	6.7	17.0	20.8	15.0	85.0	6	14.2	5.6	1.8
250153	ECX111029A	2.3	7.5		3.3	7.8	20.8		4	5.2	2.9	-7.1
250154	ECX112019A	11.0	6.8		6.2	5.0	7.0	35.9	5	7.2	2.4	-5.2
250155	ECX113019A	8.5	13.6	0.0	4.5	11.5	18.3	56.4	6	9.4	6.5	-3.0
250156	ECX114019A	10.3	9.2	14.6	9.3	20.0	15.8	79.2	6	13.2	4.4	0.8
250157	ECX115019A		7.0					7.0	1	7.0	0.0	-5.4
250158	ECX116019A	6.0	11.1	8.3	12.5	10.8	16.8	65.4	6	10.9	3.7	-1.5
250159	ECX117019A	13.8	18.5	16.0	14.0	16.5	8.5	87.3	6	14.5	3.4	2.2
250160	ECX118019A	11.2	18.0	10.8	13.9	13.8	15.1	82.7	6	13.8	2.7	1.4
250161	ECX119019A	4.0	5.8	6.3	4.0	4.0	5.9	30.0	6	5.0	1.1	-2.4
250162	ECX120019A	6.9	10.8	10.9	14.8	12.8	15.5	71.7	6	12.0	3.1	-0.4
250231	ECT111019A		7.0		10.0	22.5	14.8	54.3	4	13.6	6.8	-1.2
250232	ECT111029A	18.5		20.3	13.0	15.5		67.3	4	16.8	3.2	-4.5
250233	ECT112019A	14.3	12.0	12.5	14.0	9.8	15.8	78.3	6	13.1	2.1	0.7
250234	ECT113019A	9.5	10.0	4.8	9.3	13.5	10.0	57.1	6	9.5	2.8	-2.8
250235	ECT114019A		16.0	10.5	17.0			43.5	3	14.5	3.5	2.1
250236	ECT115019A	9.3			4.5	2.7		16.5	3	5.5	3.4	-6.8
250237	ECT116019A							0.0	0	0.0	0.0	-12.4
250238	ECT117019A	7.3			17.2	15.3	13.5	53.2	4	13.3	4.3	0.9
250239	ECT118019A					14.5		14.5	1	14.5	0.0	2.1
250240	ECT119019A				10.0			10.0	1	10.0	0.0	-2.4
250241	ECT120019A				9.5	32.0		41.5	2	20.8	15.9	8.4
251336	ECT126019A				8.2			8.2	1	8.2	0.0	-4.2
251336D	FCT126019AD	7.5		3.0			5.5	16.0	3	5.3	2.3	-7.0
Control-IR		13.5	14.5	10.0	14.5	13.5	14.7	80.7	6	13.4	1.8	1.1
Control-2R		5.7	10.0	16.0		10.5		42.2	4	10.5	4.2	-1.8
Control-3R		7.0	17.0	20.0	8.0		10.5	62.5	5	12.5	5.7	0.1
Combined								185.3	15	12.4	6.5	2.6
Pooled STD		4.21										

TABLE 10c. CARROT Average Root Lengths (cm)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR11019A	14.0	16.4	15.0	10.3	16.3	22.2	94.1	6	15.7	3.9	-4.3	3.7	YES
250022D	ECR11019AD	18.8	12.6	16.4	20.1	19.0	14.0	100.9	6	16.8	3.0	-5.4	3.7	YES
250023	ECR11029A	15.6	12.2	18.4	11.0	7.4	21.0	85.6	6	14.3	5.0	-2.8	3.7	
250024	ECR112019A	7.8	4.0	13.8	9.9	8.2	16.5	60.1	6	10.0	4.5	-1.4	3.7	
250025	ECR113019A	6.4	19.6	13.5	23.9	22.1	16.9	102.4	6	17.1	6.4	-5.6	3.7	YES
250026	ECR114019A	8.0	19.6	16.2	12.3	17.0	23.3	96.3	6	16.1	5.4	-4.6	3.7	YES
250027	ECR115019A	15.0	15.3	13.9	18.2	17.0	19.3	98.6	6	16.4	2.1	-5.0	3.7	YES
250028	ECR116019A	10.5	11.8	10.5	8.2	8.8	11.2	61.0	6	10.2	1.4	-1.3	3.7	
250029	ECR117019A	5.2	4.8	13.2	9.3	10.3	8.3	50.9	6	8.5	3.2	-2.9	3.7	
250030	ECR118019A	15.4	14.3	17.3	11.0	11.4	11.7	81.1	6	13.5	2.6	-2.1	3.7	
250031	ECR119019A	15.0	16.3	18.7	10.3	10.1	12.8	83.2	6	13.9	3.4	-2.4	3.7	
250032	ECR120019A							0.0	0	0.0	0.0	-11.4	NA	YES
250152	ECX111019A		12.5	2.8	15.3	9.6	11.0	51.2	5	10.2	4.6	-1.2	3.9	
250152D	ECX111019AD	15.5	14.4	10.0	4.9	8.3	9.9	63.0	6	10.5	3.9	-0.9	3.7	
250153	ECX111029A	9.0	12.1	11.8	16.7	5.6	12.4	67.6	6	11.3	3.7	-0.2	3.7	
250154	ECX112019A	13.0			7.3		2.8	23.0	3	7.7	3.1	-3.8	4.8	
250155	ECX113019A	20.8	10.3	11.6	22.5	10.7	11.3	87.2	6	14.5	5.5	-3.1	3.7	
250156	ECX114019A	7.3	12.9	9.2	13.6	15.5	10.4	68.9	6	11.5	3.0	-0.1	3.7	
250157	ECX115019A		5.5	13.5		15.5	7.0	41.5	4	10.4	4.9	-1.0	4.3	
250158	ECX116019A	14.0	18.9	16.8	23.8	16.5	19.2	109.2	6	18.2	3.3	-6.8	3.7	YES
250159	ECX117019A	8.8	13.7	17.6	14.7	12.8	8.3	75.9	6	12.7	3.6	-1.2	3.7	
250160	ECX118019A	6.7	11.0	13.5	7.3	12.8	7.4	58.6	6	9.8	3.0	-1.7	3.7	
250161	ECX119019A	6.3	6.8	6.3	9.6	5.5	5.8	40.2	6	6.7	1.5	-4.7	3.7	YES
250162	ECX120019A	11.0	13.4	20.0	13.8	10.6	7.5	76.3	6	12.7	4.2	-1.3	3.7	
250231	ECT111019A	13.4	16.3	17.8		13.6	17.0	78.0	5	15.6	2.0	-4.2	3.9	YES
250232	ECT111029A		10.8	16.1	31.3	16.5	11.5	86.2	5	17.2	8.2	-5.8	3.9	YES
250233	ECT112019A	10.5	11.7	14.3	13.5	12.4	17.0	79.4	6	13.2	2.3	-1.8	3.7	
250234	ECT113019A	10.0	11.9	14.2	16.9	13.8	11.8	78.5	6	13.1	2.4	-1.7	3.7	
250235	ECT114019A		7.0			31.0		38.0	2	19.0	17.0	-7.6	5.7	YES
250236	ECT115019A	10.7						10.7	1	10.7	4.4	-0.7	7.9	
250237	ECT116019A							0.0	0	0.0	0.0	-11.4	NA	YES
250238	ECT117019A					3.0	3.0	1	3.0	0.0	-8.4	7.9	YES	
250239	ECT118019A							0.0	0	0.0	0.0	-11.4	NA	YES
250240	ECT119019A							0.0	0	0.0	0.0	-11.4	NA	YES
250241	ECT120019A						9.0	9.0	1	9.0	0.0	-2.4	7.9	
251336	ECT126019A				13.0			13.0	1	13.0	0.0	1.6	7.9	
251336D	ECT126019AD			0.5	15.5			16.0	2	8.0	10.6	-3.4	5.7	
Control-1R		14.1	8.5	10.2	20.0	7.3		60.0	5	12.9	5.2	0.6	3.9	
Control-2R		12.8	18.0	10.3	11.0	16.0	11.8	79.9	6	13.3	3.0	-1.9	3.7	
Control-3R		10.0		4.0			6.0	20.0	3	6.7	3.1	-6.8	4.8	
Combined								159.9	14	11.4	4.5	0.0	2.8	
Pooled STD								4.32						

010039

December 16, 2004

Southwest Research Institute

TABLE 10d. RADISH Average Root Lengths (cm)

Sample ID	Chen ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	9.9	9.6	11.3	11.5	12.4	6.3	61.0	6	10.2	2.2	-1.7	1.5	YES
250022D	ECR111019AD	12.6	11.9	11.5	11.5	11.3	11.9	70.7	6	11.8	0.5	-0.1	1.5	
250023	ECR111029A	10.9	11.3	11.7	13.7	10.0	10.4	68.0	6	11.3	1.3	-0.6	1.5	
250024	ECR112019A	9.6	8.5	8.3	8.3	8.3	9.3	52.2	6	8.7	0.6	-3.2	1.5	YES
250025	ECR113019A	11.4	14.8	13.7	12.0	11.6	14.2	77.7	6	13.0	1.5	-1.0	1.5	
250026	ECR114019A	11.8	9.1	11.1	11.6	11.9	11.0	66.6	6	11.1	1.0	-0.8	1.5	
250027	ECR115019A	11.0	11.7	11.1	11.5	15.4	11.6	72.3	6	12.1	1.7	-0.1	1.5	
250028	ECR116019A	11.5	12.9	11.4	14.9	9.4	8.0	68.1	6	11.3	2.4	-0.6	1.5	
250029	ECR117019A	6.3	6.0	7.8	6.6	5.0	5.8	37.5	6	6.2	0.9	-5.7	1.5	
250030	ECR118019A	11.8	12.2	11.8	8.8	6.5	5.9	57.0	6	9.5	2.8	-2.4	1.5	YES
250031	ECR119019A	9.6	7.6	9.6	9.8	7.6	10.3	54.5	6	9.1	1.2	-2.8	1.5	YES
250032	ECR120019A	7.1	11.3	17.3	10.0	11.8	11.0	68.8	6	11.5	3.3	-0.4	1.5	
250152	ECX111019A	14.0	14.8	10.3	12.4	9.8	61.2	5	12.2	2.2	0.3	1.6		
250152D	ECX111019AD	12.0	12.9	10.5	10.9	11.6	10.8	68.7	6	11.5	0.9	-0.5	1.5	
250153	ECX111029A	8.4	9.3	9.0	12.3	10.1	10.9	60.0	6	10.0	1.4	-1.9	1.5	YES
250154	ECX112019A	15.5	16.0				12.0	43.5	3	14.5	2.2	2.6	2.0	YES
250155	ECX113019A	8.5	11.7	7.2	8.1	9.2	8.4	53.1	6	8.8	1.5	-3.1	1.5	YES
250156	ECX114019A	10.9	13.0	10.4	10.6	11.3	10.3	66.5	6	11.1	1.0	-0.8	1.5	
250157	ECX115019A		9.3				9.3	1	9.3	0.0	-2.7	3.4		
250158	ECX116019A	10.0	10.0	5.7	11.0	8.9	10.7	56.3	6	9.4	1.9	-2.5	1.5	YES
250159	ECX117019A	11.4	10.4	8.3	9.3	8.4	9.5	57.2	6	9.5	1.2	-2.4	1.5	YES
250160	ECX118019A	5.7	10.0	10.6	10.8	8.8	7.8	53.7	6	9.0	2.0	-3.0	1.5	YES
250161	ECX119019A	7.7	10.3	9.9	7.5	11.1	4.8	51.2	6	8.5	2.3	-3.4	1.5	YES
250162	ECX120019A	9.4	9.8	11.6	8.8	7.5	11.8	58.8	6	9.8	1.6	-2.1	1.5	YES
250231	ECT111019A	9.5	12.5	12.5	10.6	11.5	12.6	69.3	6	11.5	1.3	-0.4	1.5	
250232	ECT111029A	10.1	9.6	13.5	12.1	10.8	10.5	66.6	6	11.1	1.4	-0.8	1.5	
250233	ECT112019A	10.2	9.2	10.0	11.2	8.4	12.1	61.1	6	10.2	1.3	-1.7	1.5	YES
250234	ECT113019A	8.0	10.0	10.6	9.9	10.1	10.7	59.3	6	9.9	1.0	-2.0	1.5	YES
250235	ECT114019A	10.0	9.5	10.1	8.9	10.1	10.7	59.3	6	9.9	0.6	-2.0	1.5	YES
250236	ECT115019A	9.9	14.2	12.7	9.9	8.8	10.0	65.4	6	10.9	2.0	-1.0	1.5	
250237	ECT116019A	12.6	10.9	10.2	11.0	9.1	12.5	66.2	6	11.0	1.4	-0.9	1.5	
250238	ECT117019A	10.6	10.5	9.8	10.5	8.4	10.2	60.0	6	10.0	0.8	-1.9	1.5	YES
250239	ECT118019A	9.6	9.3		11.8	9.5	16.0	56.1	5	11.2	2.9	-0.7	1.6	
250240	ECT119019A	11.7		15.3	11.5	23.8	13.5	75.8	5	15.2	5.1	-3.3	1.6	YES
250241	ECT120019A	18.0	13.5	9.8	10.1	11.0	8.0	70.5	6	11.7	3.6	-0.2	1.5	
251336	ECT126019A	10.2	11.2	8.5	13.5	9.8	12.0	65.1	6	10.8	1.8	-1.1	1.5	
251336D	ECT126019AD	8.7	8.5	11.4	9.4	9.1	10.1	57.2	6	9.3	1.1	-2.4	1.5	YES
Control-1R		11.6	10.3	10.9	10.5	10.1	10.1	63.5	6	10.6	0.0	-1.3	1.5	
Control-2R		14.3	16.4	16.2	14.5	11.2	12.2	84.8	6	14.1	2.1	2.2	1.5	
Control-3R		13.6	11.2	10.3	11.8	9.3	10.2	66.1	6	11.0	1.5	-0.9	1.5	
Combined								214.4	18	11.9	2.2	0.0	1.1	

Pooled STD 1.88

010040

December 16, 2004

Southwest Research Institute

TABLE 10e. SOYBEAN Average Root Lengths (cm)

Sample ID	Client ID	Container 1	Container 2	Container 3	Container 4	Container 5	Container 6	Total	Number of Measurements	Mean	Standard Deviation	Difference from Control	Least Significant Difference @ 95%	Statistical Difference
250022	ECR111019A	11.7	18.0	15.0	27.7	25.3	17.2	114.8	6	19.1	6.2	-4.6	3.1	YES
250022D	ECR111019AD	15.0	14.7	16.7	16.0	18.8	30.7	111.8	6	18.6	6.1	-5.1	3.1	YES
250023	ECR111029A	17.0	17.0	16.2	26.7	30.3	30.3	137.5	6	22.9	6.9	-0.9	3.1	
250024	ECR112019A	8.3	11.8	13.3	11.0	12.0	13.7	70.2	6	11.7	1.9	-12.1	3.1	YES
250025	ECR113019A	14.7		20.2	20.0	17.8	16.0	88.6	5	17.7	2.4	-6.1	3.1	YES
250026	ECR114019A	12.0	15.0	19.0	21.0	18.5	21.5	107.0	6	17.8	3.7	-5.9	3.1	YES
250027	ECR115019A	12.2	16.3	13.8	13.2	13.5	24.3	93.3	6	15.6	4.5	-8.2	3.1	YES
250028	ECR116019A	13.5	10.0	11.8	14.8	13.8	16.8	81.2	6	13.5	2.3	-10.2	3.1	YES
250029	ECR117019A	7.8	8.7	8.7	7.5	7.7	15.0	55.3	6	9.2	2.9	-14.5	3.1	YES
250030	ECR118019A	14.3	15.0	18.7	16.8	15.7	16.2	97.3	6	16.2	1.5	-7.5	3.1	YES
250031	ECR119019A	14.5	15.8	17.7	19.3	11.3	14.0	92.7	6	15.4	2.8	-8.3	3.1	YES
250032	ECR120019A	12.7	14.8	18.0	20.8	13.2	22.2	101.6	6	16.9	4.0	-6.8	3.1	YES
250152	ECX111019A	16.0	14.8	16.7	18.0	15.3	14.0	94.8	6	15.8	1.4	-8.0	3.1	YES
250152D	ECX111019AD	13.8	14.0	13.2	15.5	15.5	15.5	87.5	6	14.6	1.0	-9.2	3.1	YES
250153	ECX111029A	13.3	14.2	14.0	16.7	13.7	15.8	87.7	6	14.6	1.3	-9.2	3.1	YES
250154	ECX112019A	10.3	13.2	13.2	12.0	9.3	5.2	63.2	6	10.5	3.0	-13.2	3.1	YES
250155	ECX113019A	13.5	14.2	19.3	16.8	15.3	16.5	95.7	6	15.9	2.1	-2.8	3.1	YES
250156	ECX114019A	16.5	16.0	16.0	14.0	18.8	13.0	94.3	6	15.7	2.9	-8.0	3.1	YES
250157	ECX115019A	13.8	21.5	24.0	6.0	13.3	19.3	97.8	6	16.3	6.6	-7.5	3.1	YES
250158	ECX116019A	10.7	16.2	15.7	13.7	15.2	15.0	86.3	6	14.4	2.0	-9.4	3.1	YES
250159	ECX117019A	10.2	17.2	24.0	19.8	23.0	26.3	120.5	6	20.1	5.8	-3.7	3.1	YES
250160	ECX118019A	7.5	7.2	14.8	10.7	10.7	11.7	62.4	6	10.4	2.8	-13.4	3.1	YES
250161	ECX119019A	10.5	13.3	11.3	14.5	20.5	12.0	82.2	6	13.7	3.6	-10.1	3.1	YES
250162	ECX120019A	12.7	15.7	14.3	15.3	16.0	89.7	89.7	6	14.9	1.3	-8.8	3.1	YES
250231	ECT111019A	12.3	15.3	13.7	13.8	14.3	18.0	87.5	6	14.6	1.9	-9.2	3.1	YES
250232	ECT111029A	11.3	12.0	14.3	13.3	15.2	12.7	78.8	6	13.1	1.4	-10.6	3.1	YES
250233	ECT112019A	14.7	12.8	17.2	13.7	15.2	17.7	91.2	6	15.2	1.9	-8.6	3.1	YES
250234	ECT113019A	14.7	14.7	15.8	15.8	14.8	19.7	95.5	6	15.9	1.9	-7.9	3.1	YES
250235	ECT114019A	14.0	9.3	7.8	8.5		11.2	50.8	5	10.2	2.5	-13.6	3.1	YES
250236	ECT115019A	12.8	18.4	23.7	15.2	13.5		83.6	5	16.7	4.5	-7.0	3.1	YES
250237	ECT116019A	9.6		15.0	12.0		13.8	50.3	4	12.6	2.3	-11.2	3.6	YES
250238	ECT117019A	13.5	12.5	13.3	14.2	14.3	15.0	82.8	6	13.8	0.9	-10.0	3.1	YES
250239	ECT118019A		36.0					36.0	1	36.0	0.0	-12.2	6.8	YES
250240	ECT119019A							0.0	0	0.0	0.0	-23.8	NA	YES
250241	ECT120019A	11.8	14.0	15.0	14.3			55.2	4	13.8	1.4	-10.0	3.6	YES
251336	ECT126019A	10.8	13.8	10.8	13.1	14.0	13.6	76.0	6	12.7	1.5	-11.1	3.1	YES
251336D	ECT126019AD	15.3	10.4	11.8	13.6	11.0	32.4	94.4	6	15.7	8.4	-8.0	3.1	YES
Control-1R		23.3	20.7	22.3	48.5	21.7	20.0	126.5	6	21.1	1.7	-2.7	3.1	
Control-2R		18.0	16.7	23.3	26.7	27.0	30.8	142.5	6	23.8	5.5	0.0	3.1	
Control-3R		18.3	26.7		35.3	31.0		111.3	4	27.8	7.3	4.1	3.6	
Combined								380.3	16	23.8	5.4	0.0	2.2	
Pooled STD		3.75												

**010041**

**Table 11. Sample Planting and Harvesting**

Sample ID	Client ID	Date Planted	Date Harvested
250022	ECR111019A	09/11/04	10/11/04
250022D	ECR111019AD	09/11/04	10/11/04
250023	ECR111029A	09/11/04	10/11/04
250024	ECR112019A	09/11/04	10/11/04
250025	ECR113019A	09/11/04	10/11/04
250026	ECR114019A	09/11/04	10/11/04
250027	ECR115019A	09/11/04	10/11/04
250028	ECR116019A	09/11/04	10/11/04
250029	ECR117019A	09/11/04	10/11/04
250030	ECR118019A	09/11/04	10/11/04
250031	ECR119019A	09/11/04	10/11/04
250032	ECR120019A	09/11/04	10/11/04
250152	ECX111019A	09/14/04	10/14/04
250152D	ECX111019AD	09/14/04	10/14/04
250153	ECX111029A	09/14/04	10/14/04
250154	ECX112019A	09/14/04	10/14/04
250155	ECX113019A	09/14/04	10/14/04
250156	ECX114019A	09/14/04	10/14/04
250157	ECX115019A	09/14/04	10/14/04
250158	ECX116019A	09/14/04	10/14/04
250159	ECX117019A	09/14/04	10/14/04
250160	ECX118019A	09/14/04	10/14/04
250161	ECX119019A	09/14/04	10/14/04
250162	ECX120019A	09/14/04	10/14/04
250231	ECT111019A	09/11/04	10/11/04
250232	ECT111029A	09/11/04	10/11/04
250233	ECT112019A	09/14/04	10/14/04
250234	ECT113019A	09/14/04	10/14/04
250235	ECT114019A	09/21/04	10/21/04
250236	ECT115019A	09/21/04	10/21/04
250237	ECT116019A	09/21/04	10/21/04
250238	ECT117019A	09/21/04	10/21/04
250239	ECT118019A	09/21/04	10/21/04
250240	ECT119019A	09/21/04	10/21/04
250241	ECT120019A	09/21/04	10/21/04
251336	ECT126019A	09/21/04	10/21/04
251336D	ECT126019AD	09/21/04	10/21/04
Control-1R		10/06/04	11/05/04
Control-2R		10/06/04	11/05/04
Control-3R		10/06/04	11/05/04

**010042**

## Appendix A

Images 1 - 6

010043

*Southwest Research Institute*

*December 16, 2004*



**Image 1.** 12 Large metal Halide HID Lamps were required for lighting (09-21-04)




**Image 2.** Tray Configuration

A-1


010044

*Southwest Research Institute*

*December 16, 2004*




Corn



Carrot

Radish



Bean

Images 3a-e. Seeds of the five test species.

**010045**

*Southwest Research Institute*

*December 16, 2004*



**Image 4a-b: Plants near midpoint (day 13) of the 30-day test (09/13/04)**

**010046**

*December 16, 2004*

*Southwest Research Institute*



Image 5a-b: Close up of plants day 13, a) healthy, b) wilting

A-4

010047

*Southwest Research Institute*

*December 16, 2004*



Image 6a-b: Controls on day 29, a) view of sample tray, b) surface of soil with green fungus

A-5

**010048**

## Appendix B

### Seed Packets

**010049**

*Southwest Research Institute*

*December 16, 2004*



**Roy's Calais Corn, Flint** 2390

True Vermont heirloom which was harvested in 1890 and 1910. Used in the Civil War to feed horses. Certified organic. Order for more Corn, certified organic.

Seeds: 99% VVOC. Lot # BPI-05-1

Corn - 1.000g = 1650 seeds, sown = 800 feet



**Seed Packet – Corn**

B-1

010050

*Southwest Research Institute*

*December 16, 2004*



**Walla Walla Sweet Onion** 2680

Famous for its sweet flavor, this is a superior variety from the Pacific Northwest. Best for eating fresh.

Germ. 89% 11/2003 Lot # TIER9-19  
B size = 1/8 ounce, ~ 728 seeds, sows ~ 50 feet



**ONIONS**




Seed Packet – Onions

B-2

010051

*Southwest Research Institute*

*December 16, 2004*




Seed Packet – Carrots

010052

*Southwest Research Institute*

*December 16, 2004*




Seed Packet - Radishes

010053

*Southwest Research Institute*

*December 16, 2004*



Seed Packet - Soybean

B-5

**010054**

## **Appendix C**

### **Chain of Custody**

December 16, 2004

435.20  
9/11/2002  
Rev. 02

## INEEL CHAIN OF CUSTODY FORM

21910

Page 1 of 1

See Instructions On Back

<sup>1</sup> Sampler (Printed): <i>Tom Haney</i>	<sup>2</sup> Sampler (Signature): <i>Tom H</i>	<sup>3</sup> Project Name: <u>Long-term Ecological Monitoring</u>						
<sup>4</sup> Laboratory Shipped To: <u>SWRI</u>		<sup>5</sup> Sampling & Analysis Plan Number: <u>ICP EXTOF - 00165</u>						
		<sup>6</sup> TOS/SOW/PSR Number: <u>473</u>						
<sup>7</sup> Sample ID#	<sup>8</sup> Sample Date	<sup>9</sup> Sample Time	<sup>10</sup> Sample Location	<sup>11</sup> Depth	<sup>12</sup> Sample Matrix	<sup>13</sup> Analysis Type No(s)	<sup>14</sup> Preservative	<sup>15</sup> Remarks
ECR111019A	X9AUG1355X	Plot	1	0-12	Soil	See Comments	INTACT	✓
ECR111029A	X03AUG1055X		1				INTACT	✓
ECR112019A	X03AUG1010X		2				INTACT	✓
ECR113019A	X9AUG1325X		3				INTACT	✓
ECR114019A	X9AUG1445X		4				INTACT	✓
ECR115019A	X09AUG0920X		5				INTACT	✓
ECR116019A	X09AUG1400X		6				INTACT	✓
ECR117019A	X09AUG1330X		7				INTACT	✓
ECR118019A	X09AUG1330X		8				INTACT	✓
ECR119019A	X09AUG1455X		9				INTACT	✓
ECR120019A	X09AUG1520X		10				INTACT	✓
Client: Bechtel SWRI Idaho, LLC IRR: 26410 Project #: 86386-38.00X Case #: ER-SOW-473 PTSR: 06/26/2004 Time: 08:38AM All samples received @ SWRI intact.								

<sup>16</sup> Comments: 9A = Earthworm Toxicity Test - Rye Grass Growth Test      Thermometer 027  
TEMPERATURE 22.0°C

Cooler Number(s):		<sup>17</sup> Relinquished By (Printed)	<sup>18</sup> Relinquished By (Signature)	<sup>19</sup> Date	<sup>20</sup> Time	<sup>21</sup> Received By (Printed)	<sup>22</sup> Received By (Signature)	<sup>23</sup> Date	<sup>24</sup> Time
		<i>Tom Haney</i>	<i>Tom H</i>	12/16/04	0900				
						<i>Norma A. Salcedo</i>	<i>Get off my back by 0830am</i>		

Distribution: Original &amp; Yellow: Accompany Shipment To Laboratory

Pink: Forward To Sample Management

Green: Retained By Project

010055

December 16, 2004

435.20  
9/11/2002  
Rev. 02

## INEEL CHAIN OF CUSTODY FORM

21911

Page 1 of 1

See Instructions On Back

<sup>1</sup> Sampler (Printed): <i>Tom Haney</i>	<sup>2</sup> Sampler (Signature): <i>[Signature]</i>	<sup>3</sup> Project Name: <i>Long-Term Ecological Monitoring</i>						
<sup>4</sup> Laboratory Shipped To:		<sup>5</sup> Sampling & Analysis Plan Number: <i>ICP-EXTOP-00165</i>						
		<sup>6</sup> TOS/SOW/PSR Number: <i>473</i>						
<sup>7</sup> Sample ID#	<sup>8</sup> Sample Date	<sup>9</sup> Sample Time	<sup>10</sup> Sample Location	<sup>11</sup> Depth (Inch)	<sup>12</sup> Sample Matrix	<sup>13</sup> Analysis Type No(s)	<sup>14</sup> Preservative	<sup>15</sup> Remarks
ECX111019A ✓	24AU009	0950✓	Plot 1	0-12	Soil ✓	9A ✓	NA ✓	INTACT ✓
ECX111029A ✓	24AU009	0950✓		1		soil ✓		INTACT ✓
ECX112019A ✓	24AU009	1045✓		2		Communt ✓		INTACT ✓
ECX113019A ✓	24AU009	1410✓		3				INTACT ✓
ECX114019A ✓	24AU009	1305✓		4				INTACT ✓
ECX115019A ✓	24AU009	1435✓		5				INTACT ✓
ECX116019A ✓	24AU009	1455✓		6				INTACT ✓
ECX117019A ✓	24AU009	1530✓		7				INTACT ✓
ECX118019A ✓	24AU009	1310✓		8				INTACT ✓
ECX119019A ✓	24AU009	1210✓		9				INTACT ✓
ECX120019A ✓	24AU009	1135✓		10				INTACT ✓

\*Comments: 9A = Earthworm Toxicity Test + Rye Grass Growth Test

Cooler Number(s):					
<sup>15</sup> Relinquished By (Printed): <i>Tom Haney</i>	<sup>16</sup> Relinquished By (Signature): <i>[Signature]</i>	<sup>17</sup> Date: <i>25AU009 0900</i>	<sup>18</sup> Time:	<sup>19</sup> Received By (Print)	
Client: Bechtel BWXT Idaho LLC SRR: 26431 Project #: 06356.38.00X Case #: ER-SOW-473 VTSR: 08/26/2004 Time: 08:30PM					

Distribution: Original &amp; Yellow: Accompany Shipment To Laboratory

Pink: Forward To Sample Management

Green: Retained By Project

*Guaranteed delivery 0830 AM*  
*THE TEMPERATURE IS 22.0°C*

010056

December 16, 2004

435.20  
9.11.2002  
Rev. 02

## INEL CHAIN OF CUSTODY FORM

21387

Page 1 of 1

See Instructions On Back

<sup>1</sup> Sampler (Printed): <b>Tom Haney</b>	<sup>2</sup> Sampler (Signature): <i>[Signature]</i>	<sup>3</sup> Project Name: <b>Long-Term Ecological Monitoring</b>						
<sup>4</sup> Laboratory Shipped To: <b>SWRI</b>		<sup>5</sup> Sampling & Analysis Plan Number: <b>ICP-EXT-07-00165</b>						
		<sup>6</sup> TOS/SOW/PSR Number: <b>473</b>						
<sup>7</sup> Sample ID#	<sup>8</sup> Sample Date	<sup>9</sup> Sample Time	<sup>10</sup> Sample Location	<sup>11</sup> Depth Inches	<sup>12</sup> Sample Matrix	<sup>13</sup> Analysis Type No(s)	<sup>14</sup> Preservative	<sup>15</sup> Remarks
ECT111019A			Plot 1	0-12	Soil	9A	None	INTACT
ECT111029A				1		Earthworm		INTACT
ECT112019A	25AUG04	10:20		2		Toxicity		INTACT
ECT113019A	25AUG04	10:45		3		Test		INTACT
ECT114019A	25AUG04	11:15		4				INTACT
ECT115019A	25AUG04	0843		5		RyeGrass		INTACT
ECT116019A	25AUG04	0915		6		Grain		INTACT
ECT117019A	25AUG04	0935		7		Test		INTACT
ECT118019A	25AUG04	0958		8				INTACT
ECT119019A	25AUG04	1025		9				INTACT
ECT120019A	25AUG04	1055		10				INTACT

<sup>16</sup> Client: Battelle BWXT Idaho, LLC  
SSRN: 264238  
Project #: 06386-38.00X  
Case #: ER-SOM-073  
VTSR: 06/27/2004 Time: 08:16AM  
All samples received @ SWRI intact.

*THERMOMETER 027  
TEMPERATURE 22.0 °C*

<sup>17</sup> Relinquished By (Printed): <b>Tom Haney</b>	<sup>18</sup> Relinquished By (Signature): <i>[Signature]</i>	<sup>19</sup> Date: <b>12/04</b>	<sup>20</sup> Time: <b>1200</b>	<sup>21</sup> Received By (Printed): <b>Garrison S. Sauerhoff</b>	<sup>22</sup> Received By (Signature): <i>[Signature]</i>	<sup>23</sup> Date: <b>06/27/04</b>	<sup>24</sup> Time: <b>08:15AM</b>
--	--	-------------------------------------	------------------------------------	--	--	--	---------------------------------------

Distribution: Original &amp; Yellow: Accompany Shipment To Laboratory

Pink: Forward To Sample Management

Green: Retained By Project

010057

December 16, 2004

435.20  
9/11/2002  
Rev 02

## INEEL CHAIN OF CUSTODY FORM

25909

Page 1 of 1

See Instructions On Back

<sup>1</sup> Sampler (Printed) <i>Tom Haney</i>	<sup>12</sup> Sampler (Signature) <i>Tom Haney</i>	<sup>1</sup> Project Name: Long-Term Ecological Monitoring	
<sup>4</sup> Laboratory Shipped To: SWRI		<sup>3</sup> Sampling & Analysis Plan Number: ICP EXT 04 - 00165	
<sup>5</sup> Sample ID# ECT 126019A	<sup>6</sup> Sample Date 8SEP04/1940	<sup>7</sup> Sample Time TSF - 07 10-12	<sup>8</sup> Sample Location Soil
<sup>9</sup> Depth inch	<sup>10</sup> Sample Matrix tox	<sup>11</sup> Analysis Type (N(s)) analyze	<sup>13</sup> Preservative Rye grass
<sup>14</sup> Remarks Earthworm NA - Intact ASAP			
<sup>15</sup> Comments			

Cust: Battelle BWXT Idaho, LLC  
SPR: 29620  
Project #: 86344-03-00X  
Case #: ER-SOW-473  
VTSR: 08/16/2004 Time: 08:20AM  
All samples received at SwRI intact.

Thermometer 027  
Temperature 22.0°C

Cooler Number(s):		<sup>12</sup> Date	<sup>13</sup> Time	<sup>14</sup> Received By (Printed)	<sup>15</sup> Received By (Signature)	<sup>16</sup> Date	<sup>17</sup> Time
<sup>11</sup> Relinquished By (Printed)	<sup>12</sup> Relinquished By (Signature)						
- Tom Haney	<i>OF-L6</i>	9SEP04	0900				

Distribution: Original &amp; Yellow: Accompany Shipment To Laboratory

Pink: Forward To Sample Management

Green: Retained By Project

010058

